

D

Delasco
Spring 2025





Website

delasco.com
Order Securely
Online 24/7



Call

800-831-6273
Talk to an Expert
M-F 8AM-5PM CST



Email

questions@delasco.com
Email Orders and
Questions to an Expert

Prices are subject to change. Visit delasco.com for the latest prices.

Contents

01 Microneedling

Microneedling Pen Kit

02-07 Chemicals

Peeling Agents, Hemostatic Agents, Mycology,
Screening & Stains

08-10 Skincare

Cleansers, Exfoliants, Anti-Aging & Moisturizers

11-21 Surgery

Punches, Curettes, Sutures, Blades, Scalpels,
Closure, Compression, Marking

22-23 Dermoscopy

DermLite® Dermatoscopes & Lumio® Viewers

24-35 Instruments

Instrument Care, Autoclaveable & Disposable

36-40 Sharps

Cannulas, Needles, Disposal & Syringes

41-44 Electrosurgery

Units, Forceps, Electrodes, Sheaths & Cautery

45-48 Cryosurgery

Units & Accessories

49-50 Examination

Accessories & Stands

51 Sterilization

Autoclaves & Ultrasonic Cleaners

52-54 Wound Care

Bandages, Wraps & Dressings

55-56 Lighting

57-61 Cosmetics

Supplies & Ear Piercing

62-65 Protection

Gloves, Laser & Non-Laser Eyewear

Symbols

Look for these symbols throughout the catalog to find new products and trusted favorites as well as special shipping, handling and availability information.



New Product

Exciting new items
in our product line

- ✔ Trusted Product: Products most often trusted by dermatologists
- ◇ Hazmat: Items greater than 1 oz in volume incur Hazmat charges by air
- Ⓜ Ground Shipping Only: May not ship by air
- ❄ Refrigerate: Refrigeration required
- ✖ No Canada: Not available in Canada
- ⊖ Non-Stock: Non-stock items are subject to longer lead times

Delasco Find our full selection of products at delasco.com



**Introducing
Delasco FDA-Cleared
Microneedling Pen Kit
The Latest Microneedling Technology
at a Fraction of the Cost**

When it comes to skin treatments, precision is everything. It's the difference between good results and exceptional ones. That's why we're proud to introduce the Delasco Microneedling Pen—an advanced, FDA-cleared device that sets a new standard in microneedling technology.

- Depth range of 0-2 mm
- 6 speed settings
- 10 needle cartridges and protective covers included in the kit
- Range of 6300-7700 RPM
- Cordless and rechargeable
- Storage base and power adapter with cable included
- Rubber grip for comfort and reduced slip
- Soft, quiet motor is less intimidating to patients



MICROPEN-NC

Item	Description	Price
MICROPEN	Microneedling Pen Kit	\$949.00
MICROPEN-NC	Needle Cartridges, 14 pin, 10/Box	\$149.99
MICROPEN-PS	Protective Sleeves, 50/Box	\$49.99

Delasco Preferred Program

Get exclusive perks and membership benefits that **save your practice time and money.**

FREE Care & Maintenance on All Delasco Surgical Instruments

FREE Shipping on Orders Over \$250


Discounted Pricing on Products & Supplies

Dedicated Support from an Expert Account Manager



**Ready to save?
Join now.**
delasco.com/preferred

Have questions? 800.831.6273 | questions@delasco.com



Chemicals – Peeling Agents



Trichloroacetic Acid

Trichloroacetic Acid solutions can be applied to penetrate the outer and middle layers of skin to remove damaged skin cells. Delasco offers the broadest range of TCA concentrations in the industry. Choose the exact strength you need at the right size for your practice.

- High purity ACS reagent grade ingredients
- Standardized weight/volume % potency

Item	Strength	(.5) oz	(1) oz	(2) oz	(PT) Pint
TCA/10/(size)	10%	-	\$23.87	\$35.24	\$71.61 G
TCA/15/(size)	15%	-	\$26.15	\$38.65	\$86.38 G
TCA/20/(size)	20%	-	\$27.28	\$42.05	\$120.43 G
TCA/30/(size)	30%	-	\$35.24	\$53.42	\$136.94 G
✓ TCA/35/(size)	35%	-	\$38.65	\$57.97	\$134.13 G
TCA/40/(size)	40%	-	\$39.78	\$57.97	\$143.22 G
TCA/50/(size)	50%	-	\$43.20	\$65.93	-
TCA/60/(size)	60%	-	\$45.47	\$69.33	-
TCA/75/(size)	75%	-	\$46.60	\$73.88	-
TCA/80/(size)	80%	NEW	\$42.66	\$50.01	\$76.16
TCA/85/(size)	85%	-	\$50.01	\$76.16	-
TCA/90/(size)	90%	-	\$50.01	\$76.16	-
TCA/100/(size)	100%	-	\$50.01	\$76.16	-



Dichloroacetic Acid

- Also known as bichloroacetic acid
- High purity reagent grade ingredients, 100% DCA

Item	Size	Price
DCA/1	1oz	\$62.52
DCA/2	2oz	\$89.80 ◇



Monochloroacetic Acid

- High purity ACS reagent grade ingredients
- Saturated solution
- 79% w/w = 104% w/v

Item	Size	Price
MCA/SAT/1	1oz	\$77.29
MCA/SAT/2	2oz	\$105.70 ◇



Acetone

- USP grade

Item	Size	Price
ACETONE	2oz	\$6.83 G
ACETONE/8	8oz	\$11.41 G



Jessners Solution

Formulated with the original Jessners formula containing 14% salicylic acid, 14% resorcinol, and 14% lactic acid, in an ethanol base.

- USP grade ingredients

Item	Size	Price
JS/1	1oz	\$30.69
JS	2oz	\$47.74 G
✓ JS/4	4oz	\$75.02 G
JS/PT	PT	\$213.69 G

Modified Jessners Solution

A variation on the original Jessners formula. Modified Jessners is prepared with 17% salicylic acid, 17% lactic acid, and 8% citric acid in an alcohol solution.

- No resorcinol
- USP grade ingredients

Item	Size	Price
MJS/1	1oz	\$30.69
MJS	2oz	\$47.74 ◇
MJS/4	4oz	\$75.02 ◇
MJS/PT	PT	\$213.69 ◇



Lactic Acid

- USP Grade
- 88-92%, w/w

Item	Size	Price
LA/4	4oz	\$27.28 ◇



Mandelic Acid

- 30% in Ethanol

Item	Size	Price
MAND/30/1	1oz	\$43.20 G



Phenol, Liquefied

- Minimum 89% by weight
- USP grade

Item	Size	Price
PHENOL/1	1oz	\$75.17 G
PHENOL/2	2oz	\$121.24 G

Chemicals – Peeling Agents

Glycolic Acid

Glycolic Acid is the smallest molecule of the alpha hydroxy acids. Solutions are prepared in standardized 35-70% w/v aqueous formulations.

- Cosmetic grade
- Un-buffered



Neutralizing Spray

- Sodium Bicarbonate 10%
- Neutralizes glycolic acid
- USP grade ingredients

Item	Strength	(1) oz	(2) oz	(4) oz	(8) oz
G35(size)	35%	\$43.20	\$57.97 ◇	\$86.38 G	\$129.58 G
G50(size)	50%	\$55.69	\$77.29 ◇	\$113.66 G	\$165.95 G
<input checked="" type="checkbox"/> G70(size)	70%	\$62.52	\$93.21 ◇	\$148.90 G	\$225.06 G

Item	Size	Price
GLY/N	8oz	\$17.05



Peel Prep

- Dissolves and removes lipids and oils

Item	Size	Price
<input checked="" type="checkbox"/> PP/1	1oz	\$7.96
<input checked="" type="checkbox"/> PP/4	4oz	\$12.51 G
<input checked="" type="checkbox"/> PP/PT	PT	\$25.01 G

Salicylic Acid

This popular beta hydroxy acid is available in 10-30% concentrations in a hydro-alcoholic solution and full strength in polyethylene glycol.

- USP grade ingredients

Item	Size	Price
10% with Hydroethanolic Solution		
SAHA/10/1	1oz	\$24.15
SAHA/10	2oz	\$36.37 ◇
SAHA/10/PT	PT	\$164.81 ◇
20% with Hydroethanolic Solution		
SAHA/20/1	1oz	\$24.15
SAHA/20	2oz	\$36.37 ◇
SAHA/20/PT	PT	\$164.81 ◇
30% with Hydroethanolic Solution		
SAHA/1	1oz	\$27.28
<input checked="" type="checkbox"/> SAHA	2oz	\$42.05 ◇
SAHA/PT	PT	\$211.42 ◇

Also available in 30% with Polyethylene Glycol.



Croton Oil

- Tigilium Oil

Item	Size	Price
<input checked="" type="checkbox"/> CROT/1	8ml	\$133.28
<input checked="" type="checkbox"/> CROT/2	30ml (1oz)	\$390.88



Lacti-Blend Acid Peel

- 35% Lactic Acid
- 7% Salicylic Acid
- Evens pigmentation
- Excellent for acid rotation
- Must be neutralized

Item	Size	Price
SC-120	4oz	\$60.50

Glyco-Blend Acid Peel

- 35% Glycolic Acid
- 7% Salicylic Acid
- Evens pigmentation
- Excellent for acid rotation
- Must be neutralized

Item	Size	Price
SC-121	4oz	\$60.50

Glyco-Bright Acid Peel

- 45% Glycolic Acid
- 2% Hydroquinone
- Reduces Hyperpigmentation
- Must be neutralized

Item	Size	Price
SC-123	4oz	\$60.50

Lite Jessner's Acid Peel

- Resorcinol, Lactic and Salicylic Acids are 7% by volume
- Significant exfoliation, yet slightly less penetrating than the standard Jessner's Peel
- Self neutralizing

Item	Size	Price
SC-117	4oz	\$60.50

Ultra Jessner's Acid Peel

- Standard Jessner's formula plus 2% Hydroquinone
- Self neutralizing

Item	Size	Price
SC-119MD	4oz	\$73.70

Lacti-Bright Acid Peel

- 45% Lactic Acid
- 2% Hydroquinone
- Reduces Hyperpigmentation
- Must be neutralized

Item	Size	Price
SC-122	4oz	\$60.50

Lite Modified Jessner Acid Peel

- Does not contain resorcinol
- Lactic, Salicylic, and Glycolic Acids are 7% by volume
- Self neutralizing

Item	Size	Price
SC-115	4oz	\$60.50

Lite Jessner's Bright Acid Peel

- 2% Hydroquinone
- Standard Light Jessner's formula plus 2% Hydroquinone
- Self neutralizing

Item	Size	Price
SC-117B	4oz	\$64.90

Chemicals – Hemostatic Agents



Aluminum Chloride

Aluminum Chloride solutions prepared from USP grade materials. Formulated in ethanol, isopropanol or in aqueous solutions.
 • USP grade ingredients

Item	(1) oz	(2) oz	(4) oz	(8) oz	(PT) Pint
✓ AC/20/(size) 20% in ethyl alcohol	\$20.30	\$31.13 ◇	\$42.89 ^{NEW}	\$81.84 ^{NEW} ◇	\$155.88 ◇
✓ AC/35/(size) 35% in purified water	\$13.64	\$20.46 ◇	\$34.21 ^{NEW} ◇	-	\$113.66 ◇
AC/35/70/(size) 35% in 70% isopropanol	\$23.87	\$31.83 ◇	-	-	\$179.59 ◇
AC/50/(size) 50% in 50% isopropanol	\$19.32	\$30.69 ◇	-	-	\$185.27 ◇
AC/70/(size) 70% in purified water	\$19.32	\$30.69 ◇	-	-	\$185.27 ◇



Ferric Chloride

- High purity ACS reagent grade ingredients
- 3.8 molal 37.5% w/v
- 3.8 Molar 61.6% w/v
- Reduced tattooing risk compared to Monsel's

Item	Size	Price
FC/1/1	3.8m, 1oz	\$13.64
FC/1	3.8m, 2oz	\$19.32 ◇
FC/2/1	3.8M, 1oz	\$15.92
FC/2	3.8M, 2oz	\$23.87 ◇

Monsel's Solution

- Ferric subsulfate solution

Item	Size	Price
MONSELS/1	1oz	\$13.64
MONSELS	2oz	\$21.60
MONSELS/PT	16oz	\$70.48

Tissue Marking Dyes

- Use to mark the edges of excised tissue and bases of multiple lesions when put into the same container
- Stains survive formalin fixation and H & E staining
- Compare to Davidson stains (excludes Brilliant Green)



Item	Description	(/1) ½ oz	(/2) 1oz	(/8) 8oz
BLK/(size)	Black	\$9.09	\$11.91	\$73.83
✓ BLU/(size)	Blue	\$8.66	\$12.51	\$70.48
RED/(size)	Red	\$8.66	\$12.51	\$67.12
YEL/(size)	Yellow	-	\$11.91	\$70.48
GRN/(size)	Green	-	\$11.91	\$67.12
BG/5/1	Brilliant Green 5% in Denatured Ethyl Alcohol		\$27.28	



Chlorazol Fungal Stain

- Stains hyphae green against a grey background
- Contains KOH 5% with DMSO and Chlorazol Black E contrast dye
- May examine with bright/reduced light

Item	Size	Price
C/1	1/2oz	\$19.32
C/2	1oz	\$30.69



Lugol's Solution

- USP grade ingredients
- 10% potassium iodide and 5% iodine in water

Item	Size	Price
LUGOL/1	5% Iodine, 1oz	\$21.98
LUGOL	5% Iodine, 2oz	\$31.66



Schiller's Solution

- USP grade ingredients
- 2% potassium iodide and 1% iodine in water

Item	Size	Price
SCHILLER	4oz	\$28.06



Potassium Hydroxide Preparations

- Also known as KOH
- The addition of DMSO speeds penetration for rapid use
- DMSO Heating not required

Item	Description	(1) oz	(2) oz	(PT) Pint
KOHD/10/(size)	KOH 10% with DMSO	\$18.40	\$29.23	◆
KOHD/20/(size)	KOH 20% with DMSO	\$18.40	\$29.23	◆
<input checked="" type="checkbox"/> KOH/10/(size)	KOH 10% in Water	\$12.99	\$15.59	◆ \$80.11 ◆
KOH/20/(size)	KOH 20% in Water	\$14.07	\$19.49	◆



Acetic Acid

- USP grade ingredients

Item	Strength	(4) oz	(8) oz	(PT) Pint
AA/025/(size)	0.25%	-	-	\$23.08 NEW
AA/3/(size)	3%	\$15.01	\$18.40	\$18.37
AA/4/(size)	4%	\$14.07		\$28.15
<input checked="" type="checkbox"/> AA/5/(size)	5%	\$18.36	\$20.57	\$16.26



DTM Agar Slant, Dermatophyte Test Media

- Shatter Resistant Polycarbonate tubes
- Antibiotics suppress bacterial growth
- Turns red with growth of dermatophytes
- No color changes with growth of saprophytes

Item	Size	Price
<input checked="" type="checkbox"/> DTM20	10 ml tube with 7.5 ml media on the slant, 20/box	\$117.16



Facial Cleansers

Formulated to deeply cleanse pores and remove excess oils and impurities from the skin. The Gentle Cleansing Gel works great for sensitive and dry skin. The AHA Rejuvenating Cleanser has Alpha and Beta Hydroxy Acids that gently exfoliate dead skin cells. The Purifying Facial Wash completely removes oils, and prevents pore clogging while protecting healthy skin cells.

Item	Description	Price
ST-113	Gentle Cleansing Gel, 4oz	\$20.63
STPRO-113B	Gentle Cleansing Gel, 16oz	\$47.30
ST-114	AHA Rejuvenating Cleanser, 4oz	\$20.63
STPRO-114B	AHA Rejuvenating Cleanser, 16oz	\$47.30
ST-111	Purifying Facial Wash, 4oz	\$20.63
STPRO-111B	Purifying Facial Wash, 16oz	\$47.30



Exfoliation Treatments

These exfoliation treatments are to be used as part of a home regimen to keep dead skin cells from dulling skin and enhance the results from chemical peels. Lactic Cellular Turnover helps skin retain moisture while exfoliating. Glycolic Cellular Turnover enhances exfoliation and penetrates deeper than other acids in the skin. Anti-Oil Cellular Turnover combines lactic and salicylic acid that helps clean pores while exfoliating. Mandelic Cellular Turnover gently exfoliates and leaves skin glowing. Refining Cellular Turnover combines glycolic and salicylic acid that helps reach deep cystic-like blemishes.

Item	Description	Price
ST-115	Lactic Cell Turnover, 1.75oz	\$21.18
ST-117	Glycolic Cell Turnover, 1.75oz	\$21.18
ST-119	Anti-Oil Cell Turnover, 1.75oz	\$21.18
ST-120	Mandelic Cell Turnover, 1.75oz	\$23.10
ST-121	Refining Cell Turnover, 1.75oz	\$23.10



Cleansing Pads

360 Cleansing Pads are easy to use and are two-sided for a thorough clean. The scrubby side of the pad cleans powerfully, even tough eyeliner, without irritating the skin or eyes. And no rinsing! You don't have to skip cleansing anymore after a long day.

Item	Description	Price
ST-149	360 pads, 25 applications	\$15.40



Tri-Ex Peel Pads

A triple acid exfoliation formula in a convenient pad penetrates the dermal layer to simulate new cell growth and leave the skin smooth, tight, and luminous.

Item	Size	Price
ST-143B	4oz	\$21.45



Micro-Dermal Scrub

Micro-Dermal Scrub warms like the tropical sun on your skin while it gently exfoliates. This self-warming scrub invigorates the skin and captivates your senses with its citrus scent. Retinol and peptides infuse nourishment, while dead skin cells and contaminants are removed. For a facial regimen, this is a "take me away" kind of product that you will love.

Item	Size	Price
PC-102	4oz	\$24.75



Vita-C 10

Vita-C 10 uses a protein encapsulated form of Vitamin C that is released once it's applied to the skin. Nourishment is then provided to the intercellular space that strengthens collagen fibers. Vitamin C is crucial in reducing brown spots and sunburn while improving skin tone and texture.

**This product is light and heat sensitive. Store in a temperature-controlled location out of direct sunlight.*

Item	Size	Price
ST-099A	0.5oz	\$20.90

VitaFuse Serum

Free radicals are all around us every day causing damage to our skin, but they don't stand a chance against this thick, vitamin-rich serum. The blend of Pro-Vitamin B-5, A, E, K, and Carrot Seed Extract come together to repair damage to skin associated with free radicals.

Item	Size	Price
ST-103	0.5oz	\$26.40
STPRO-103C	8oz	\$184.80



VitaFuse Lotion

This light, vitamin-rich lotion hydrates while it restores. Pro-Vitamin B-5, along with vitamins A, K, and Carrot Seed Extract bring nourishment to the skin. OPC has also been added to this blend to reduce dark circles, bruising, redness and puffiness. It's an essential part of a post-clinical treatment and a regular skincare regimen.

Item	Size	Price
ST-102	1oz	\$35.20
STPRO-102C	8oz	\$61.60



HA-B5 Hyaluronic Moisturizer

This high efficiency hyaluronic and Pro Vitamin B5 hydrator glides on and boosts the skin's ability to hold moisture while leaving no oily residue. The prominence of fine lines and wrinkles is reduced, and it provides fullness back to the skin.

Item	Size	Price
ST-148	1oz	\$35.20
STPRO-148B	4oz	\$115.50



HydraCell Turnover Booster

This plant-derived squalane serum has been shown to increase cellular turnover by 34%. 34%! Apply HydraCell after a Skin Therapeutic Cellular Turnover Treatment and push delayed gratification out the window.

Item	Size	Price
ST-144	1oz	\$18.15
STPRO-144B	8oz	\$55.55



3D Peptide Crème

3D Peptide Crème brings three powerful peptides together in one formulation to significantly improve skin's thickness and smooth fine lines and wrinkles. Consistent use will improve the tonicity of skin and reveal a more smooth sculpted appearance. The effects of lost collagen can be noticeably recovered with this powerful tri-peptide combination.

Item	Size	Price
ST-140	1oz	\$26.40



Syner-G Repair Crème

The signs of aging won't be there long with Syner-G around. Infused within this crème is a powerful synthetic peptide formula that mimics the muscle relaxation peptide found in the venom of the Temple Viper. This peptide is so effective that in vivo studies have shown a 52% reduction in wrinkle size after 28 days of use.

Item	Size	Price
ST-142	1oz	\$30.25



Radiance Crème

Radiance Crème is a tightening complex designed to reduce the appearance of fine lines and wrinkles and restore firmness to the face, neck, and décolleté. Aurasphere and Squalane work together to penetrate and fill deep lines and wrinkles, even skin tone and hydrate skin.

Item	Size	Price
ST-145	1.75oz	\$29.43



Brightinol Dark Spot Corrector

Melasma and other brown skin discolorations resulting from acne, sun damage or hormonal changes can be difficult to reverse. Brightinol Dark Spot Corrector is a powerful formulation with Retinol and Kojic Acid designed to reduce the appearance of dark spots without irritating or drying the skin.

Item	Size	Price
PC-105	0.5oz	\$28.05





DermaBoost Retinol Therapy

Collagen and elastin in the lower layer of skin are the foundation for skins' firmness and play a vital in its ability to renew cells. DermaBoost Retinol Therapy is formulated to penetrate deep and boost collagen synthesis so skin is strengthened from the bottom layer up.

Item	Size	Price
PC-106	0.5oz	\$27.80

Tinted Moisturizer SPF 30

Zinc Oxide and Titanium Dioxide provide UVA and UVB protection that absorbs quickly and has a sheer finish with no residue! Fresh scent and non-comedogenic. Chamomile Flower Extract, Cucumber Fruit Extract and Sodium Hyaluronate add beneficial nourishing properties that soothe skin, help retain moisture.

Item	Size	Price
ST-155	1.7oz	\$30.45



Post Procedure Lotion

Made with a lightweight formulation designed to provide long lasting hydration for all skin types. It is ideally used following any facial procedure or laser treatment to promote rapid recovery and healing. Acts as a moisture shield and is a perfect every day moisturizer that absorbs rapidly leaving the skin supple and smooth with no oily residue.

Item	Size	Price
ST-135	8oz	\$44.00
STPRO-135B	16oz	\$93.50



Intensive Eye Crème

The delicate skin around the eyes needs intense care. Intensive Eye Crème brings a host of essential vitamins, nutrients, proteins, and antioxidants together to create a rich hydration barrier around the eyes. Shea Butter and Hemp Seed Oil blend together to aid in diminishing the fine lines around eyes where crows feet form. This formulation also minimizes dark circles and promotes elasticity and firmness.

Item	Size	Price
ST-104B	0.5oz	\$33.55



iLift-R Treatment

iLIFT-R Treatment is a blend of Hyaluronic Acid, natural fruit extracts and vegetable proteins that work together to strengthen the delicate skin under the eyes, reduce puffiness and the appearance of dark circles while diminishing crow's feet.

Item	Size	Price
ST-147	0.5oz	\$30.80



Tri-Ex Spot Treatment

This unique blend of Salicylic, Glycolic and Lactic acids can be placed directly on a blemish with the small roller ball applicator. Apply at the first sign of a blemish and keep these unwanted visitors behind closed doors.

Item	Size	Price
PC-104	0.35oz	\$14.03

Delasco Biopsy Punches

Superior quality at an economical price!

- Razor sharp seamless stainless steel blade with a safety cap
- Unique hollow core hex-handle improves tactile feel, yet minimizes weight
- Clearly marked with punch size
- Individually packaged for single use
- Sterile and disposable

Delasco Biopsy Punch Guarantee: Try a box of Delasco Biopsy Punches. If you aren't completely satisfied, call Delasco for full credit.

Item	Description	Price
DBP-15	1.5 mm, 50/box	\$80.25
DBP-2	2 mm, 50/box	\$80.25
DBP-25	2.5 mm, 50/box	\$80.25
DBP-3	3 mm, 50/box	\$80.25
DBP-35	3.5 mm, 50/box	\$80.25
✓ DBP-4	4 mm, 50/box	\$80.25
DBP-5	5 mm, 50/box	\$80.25
DBP-6	6 mm, 50/box	\$80.25
DBP-8	8 mm, 50/box	\$80.25
DBP-A	(10) 2 mm, (15) 3, 4 mm and (5) 5, 6 mm, 50/box	\$80.25



Delasco Dermal Curettes

Superior quality at an economical price!

- Sharp, seamless, stainless steel edge
- Clearly marked with size
- Individually packaged for single use
- Sterile and disposable

Delasco Dermal Curette Guarantee: Try a box of Delasco Dermal Curettes. If you aren't completely satisfied, call Delasco for full credit.

Item	Description	Price
DDC-2	2 mm, 50/box	\$84.54
DDC-3	3 mm, 50/box	\$84.54
✓ DDC-4	4 mm, 50/box	\$84.54
DDC-5	5 mm, 50/box	\$84.54
DDC-7	7 mm, 50/box	\$84.54




Miltex® Disposable Punches

- Seamless, stainless steel, razor-sharp cutting edge with ribbed handle
- Sterile and disposable
- CE marked

Sizes available: 1, 2, 2.5, 3, 3.5, 4, 5, 6 and 8 mm.


Item	Description	Price
PUNCH/M/(size)	50/box	\$151.98




Miltex® Dermal Curettes

- Seamless, stainless steel, semi-sharp cutting edge
- Symmetrically ground beveled edge
- Clearly marked with size
- Individually packaged
- Sterile and disposable
- CE marked

Item	Description	Price
33-52	2 mm, 50/box	\$164.02
33-53	3 mm, 50/box	\$164.02
33-54	4 mm, 50/box	\$164.02
33-55	5 mm, 50/box	\$164.02
33-57	7 mm, 50/box	\$164.02





NEW

DemeTech Sutures

ABSORBABLE

Surgical Plain Gut

18"

Item	Make	Needle Style	Needle Size	Suture Size	Price
1916D	DEMETECH	PC-1	13 mm	6-0	\$101.98
1915D	DEMETECH	DPC-1	13 mm	5-0	\$98.67
686D	DEMETECH	P-3	13 mm	5-0	\$92.21
1644D	DEMETECH	P-3	13 mm	4-0	\$98.67

30"

550D	DEMETECH	DFS-2	19 mm	5-0	\$59.37
551D	DEMETECH	DFS-2	19 mm	4-0	\$60.61
554D	DEMETECH	DFS-2	19 mm	3-0	\$50.77

Surgical Chromic Gut

18"

1654D	DEMETECH	DP-3	13 mm	4-0	\$68.34
687D	DEMETECH	DP-3	13 mm	5-0	\$72.43
1635D	DEMETECH	DPS-3	16 mm	6-0	\$101.98

30"

559D	DEMETECH	DFS-2	19 mm	4-0	\$51.63
------	----------	-------	-------	-----	---------

Vicryl Rapide - Polyglactin 910

18" undyed

FA VQ493D	DEMETECH	DP-3	13 mm	5-0	\$127.99
VQ494D	DEMETECH	DP-3	13 mm	4-0	\$131.58

PGCL Monofilament 25

18" undyed

Y493D	DEMETECH	DP-3	13 mm	5-0	\$98.69
Y494D	DEMETECH	DP-3	13 mm	4-0	\$96.12
Y495D	DEMETECH	DPS-2	19 mm	5-0	\$89.26
Y496D	DEMETECH	DPS-2	19 mm	4-0	\$95.04
Y497D	DEMETECH	DPS-2	19 mm	3-0	\$94.50

18" Violet

Y463D	DEMETECH	DP-3	13 mm	5-0	\$96.18
Y464D	DEMETECH	DP-3	13 mm	4-0	\$89.26

Vicryl - Polyglactin 910

18" Undyed Braided

490D	DEMETECH	DP-1	11 mm	5-0	\$101.01
492D	DEMETECH	DP-3	13 mm	6-0	\$90.66
493D	DEMETECH	DP-3	13 mm	5-0	\$99.29
494D	DEMETECH	DP-3	13 mm	4-0	\$99.29
495D	DEMETECH	DPS-2	19 mm	5-0	\$99.53
496D	DEMETECH	DPS-2	19 mm	4-0	\$101.31
497D	DEMETECH	DPS-2	19 mm	3-0	\$97.08

18" Violet Braided

464D	DEMETECH	DP-3	13 mm	4-0	\$98.11
392D	DEMETECH	DFS-2	19 mm	4-0	\$61.26
393D	DEMETECH	DFS-2	19 mm	3-0	\$61.26

27" Undyed Braided

J443D	DEMETECH	DFS-1	24 mm	2-0	\$61.26
-------	----------	-------	-------	-----	---------



Polydioxanone

18" Clear Mono

Item	Make	Needle Style	Needle Size	Suture Size	Price
Z494D	DEMETECH	DP-3	19 mm	3-0	\$135.10
Z497D	DEMETECH	DP-3	13 mm	4-0	\$127.99

30" Violet

M397D	DEMETECH	DFS-2	19 mm	4-0	\$62.76
M452D	DEMETECH	DFS-1	24 mm	3-0	\$65.98

NON-ABSORBABLE

Polypropylene, Blue Mono

18"

8697D	DEMETECH	DP-1	11 mm	6-0	\$105.26
8695D	DEMETECH	DP-3	13 mm	6-0	\$108.74
8698D	DEMETECH	DP-3	13 mm	5-0	\$105.26
8699D	DEMETECH	DP-3	13 mm	4-0	\$105.26
8687D	DEMETECH	DPS-2	19 mm	3-0	\$94.95
8686D	DEMETECH	DPS-2	19 mm	5-0	\$95.28
8682D	DEMETECH	DPS-2	19 mm	4-0	\$98.33
8683D	DEMETECH	DFS-2	19 mm	4-0	\$62.18
8665D	DEMETECH	DFS-2	19 mm	3-0	\$56.77
8684D	DEMETECH	DFS-1	24 mm	3-0	\$53.33

10"

1011D	DEMETECH	DC-3	13 mm	6-0	\$40.44
1013D	DEMETECH	DC-3	13 mm	5-0	\$39.84
1031D	DEMETECH	DC-3	13 mm	4-0	\$39.50

18" Silk, Black Braided

1265D	DEMETECH	DP-3	13 mm	5-0	\$92.83
685D	DEMETECH	DFS	18 mm	6-0	\$47.25

Nylon, Black Mono

18"

1698D	DEMETECH	DP-3	13 mm	6-0	\$82.35
698D	DEMETECH	DP-3	13 mm	5-0	\$74.77
699D	DEMETECH	DP-3	13 mm	4-0	\$72.62
1666D	DEMETECH	DPS-2	19 mm	5-0	\$78.78
1667D	DEMETECH	DPS-2	19 mm	4-0	\$58.72
1669D	DEMETECH	DPS-2	19 mm	3-0	\$65.68
918D	DEMETECH	DFS-2	19 mm	5-0	\$42.74
922D	DEMETECH	DFS-2	19 mm	4-0	\$38.45
928D	DEMETECH	DFS-2	19 mm	3-0	\$36.41
925D	DEMETECH	DFS-1		3-0	\$41.43

Biopsy

10" Nylon, Black Mono

913D	DEMETECH	DC-3	13 mm	5-0	\$35.23
911D	DEMETECH	DC-3	13 mm	6-0	\$39.69
1273D	DEMETECH	DP-3	13 mm	6-0	\$63.87


FA Sutures are fast absorbing.

*Sutures with item numbers ending in H are packaged in boxes of 36. All other sutures are sold in boxes of 12.


Surgery – Sutures

At Delasco, we carry high-quality suture brands and materials in a wide range of sizes to suit your specific needs. For a full list of sutures, visit delasco.com.


Suture Needle Chart




P-1 (3/8 Circle)
11 mm PP
Reverse Cutting



P-3 (3/8 Circle)
13 mm PP
Reverse Cutting



PS-2 (3/8 Circle)
19 mm PP
Reverse Cutting



C-2 (3/8 Circle)
12 mm
Reverse Cutting

NON-ABSORBABLE Polypropylene, Blue Mono

18"

Item	Make	Needle Style	Needle Size	Suture Size	Price
8696G	ETHICON	P-1	11 mm	7-0	\$390.98
8697C	CP MEDICAL	NP-1	11 mm	6-0	\$115.79
8695G	ETHICON	P-3	13 mm	6-0	\$378.65
8695C	CP MEDICAL	P-3	13 mm	6-0	\$119.61
8698G	ETHICON	P-3	13 mm	5-0	\$305.43
8698C	CP MEDICAL	NP-3	13 mm	5-0	\$121.58
8699G	ETHICON	P-3	13 mm	4-0	\$352.08
✓ 8699C	CP MEDICAL	P-3	13 mm	4-0	\$121.58
8680G	ETHICON	PS-3	16 mm	6-0	\$366.39
8687H*	ETHICON	PS-2	19 mm	3-0	\$1,038.27
8687C	CP MEDICAL	NPS-2	19 mm	3-0	\$109.67
8686G	ETHICON	PS-2	19 mm	5-0	\$285.06
8686C	CP MEDICAL	NPS-2	19 mm	5-0	\$110.05
8682G	ETHICON	PS-2	19 mm	4-0	\$345.75
8682C	CP MEDICAL	NPS-2	19 mm	4-0	\$113.57
8683G	ETHICON	FS-2	19 mm	4-0	\$173.45
8683C	CP MEDICAL	NFS-2	19 mm	4-0	\$71.82
8665G	ETHICON	FS-2	19 mm	3-0	\$173.45
8665C	CP MEDICAL	NFS-2	19 mm	3-0	\$65.57
8684C	CP MEDICAL	NFS-1	24 mm	3-0	\$61.59
⊖ 8663G	ETHICON	PS-1	24 mm	3-0	\$350.00

14"

⊖ 8571G	ETHICON	TS	48mm	2-0	\$382.78
---------	---------	----	------	-----	----------

10"

1011C	CP MEDICAL	NC-3	13 mm	6-0	\$46.70
1013C	CP MEDICAL	NC-3	13 mm	5-0	\$46.01
1031C	CP MEDICAL	NC-3	13 mm	4-0	\$45.62

18" Undyed Mono

⊖ MCP494G	ETHICON	P-3	13 mm	4-0	\$419.16
⊖ MCP496G	ETHICON	PS-2	19 mm	4-0	\$392.13



Biopsy 9" Nylon, Black Mono

✓ BDL-40B	DELASCO	FS-3	16 mm	4-0	\$37.67
✓ BDL-50B	DELASCO	FS-3	16 mm	5-0	\$37.67
✓ BDL-60B	DELASCO	FS-3	16 mm	6-0	\$37.67
✓ MS-40B	DELASCO	C-3	13 mm	4-0	\$37.67
✓ MS-50B	DELASCO	C-3	13 mm	5-0	\$37.67
✓ MS-60B	DELASCO	C-3	13 mm	6-0	\$37.67

Delasco Biopsy Sutures (shown at bottom of pg. 14) are smaller than typical sutures. This minimizes tissue trauma and scarring, while providing the same flexibility and strength. Which is especially important when performing biopsies on visible areas like the face.

Suture Needle Chart

			
FS (3/8 Circle) 26 mm Reverse Cutting	FS-1 (3/8 Circle) 24 mm Reverse Cutting	FS-2 (3/8 Circle) 19 mm Reverse Cutting	FS-3 (3/8 Circle) 16 mm Reverse Cutting

18" Silk, Black Braided

Item	Make	Needle Style	Needle Size	Suture Size	Price
1265C	CP MEDICAL	NP-3	13 mm	5-0	\$107.22
736G	ETHICON	C-3	13 mm	4-0	\$121.17

Nylon, Black Mono

18"					
Item	Make	Needle Style	Needle Size	Suture Size	Price
1698G	ETHICON	P-3	13 mm	6-0	\$230.89
1698C	CP MEDICAL	NP-3	13 mm	6-0	\$95.12
698G	ETHICON	P-3	13 mm	5-0	\$230.49
698C	CP MEDICAL	P-3	13 mm	5-0	\$86.36
699G	ETHICON	P-3	13 mm	4-0	\$258.29
699C	CP MEDICAL	P-3	13 mm	4-0	\$83.87
1666C	CP MEDICAL	NPS-2	19 mm	5-0	\$90.99
1667G	ETHICON	PS-2	19 mm	4-0	\$245.84
1667C	CP MEDICAL	NPS-2	19 mm	4-0	\$67.82
1669H*	ETHICON	PS-2	19 mm	3-0	\$732.85
1669C	CP MEDICAL	PS-2	19 mm	3-0	\$75.86
661G	ETHICON	FS-2	19 mm	5-0	\$128.16
918C	CP MEDICAL	NFS-2	19 mm	5-0	\$49.36
662G	ETHICON	FS-2	19 mm	4-0	\$116.98
922C	CP MEDICAL	NFS-2	19 mm	4-0	\$44.42
928C	CP MEDICAL	NFS-2	19 mm	3-0	\$42.05
663G	ETHICON	FS-1	24 mm	3-0	\$122.16
925C	CP MEDICAL	NFS-1	24 mm	3-0	\$47.85

ABSORBABLE

Surgical Plain Gut

18"					
Item	Make	Needle Style	Needle Size	Suture Size	Price
1916C	CP MEDICAL	NPC-1	13 mm	6-0	\$117.79
1915G	ETHICON	PC-1	13 mm	5-0	\$428.41
1915C	CP MEDICAL	NPC-1	13 mm	5-0	\$113.97
686C	CP MEDICAL	NP-3	13 mm	5-0	\$106.50

30"

550C	CP MEDICAL	FS-2	19 mm	5-0	\$68.58
551C	CP MEDICAL	NFS-2	19 mm	4-0	\$70.00
554C	CP MEDICAL	NFS-2	19 mm	3-0	\$58.64

Surgical Chromic Gut

18"					
Item	Make	Needle Style	Needle Size	Suture Size	Price
1654C	CP MEDICAL	P-3	13 mm	4-0	\$78.93
687C	CP MEDICAL	P-3	13 mm	5-0	\$83.65

30"					
Item	Make	Needle Style	Needle Size	Suture Size	Price
559C	CP MEDICAL	FS-2	19 mm	4-0	\$59.63

FA Sutures are fast absorbing.

*Sutures with item numbers ending in H are packaged in boxes of 36. All other sutures are sold in boxes of 12.

Surgery – Sutures

ABSORBABLE

Vicryl Rapide - Polyglactin 910

18" undyed

Item	Make	Needle Style	Needle Size	Suture Size	Price
FA VR493	ETHICON	P-3	13 mm	5-0	\$425.51
FA VR494	ETHICON	P-3	13 mm	4-0	\$423.23
FA VQ494C	CP MEDICAL	NP-3	13 mm	4-0	\$151.98

PGCL Monofilament 25

18" undyed

Y492G	ETHICON	P-3	13 mm	6-0	\$392.26
Y493G	ETHICON	P-3	13 mm	5-0	\$387.53
Y493C	CP MEDICAL	NP-3	13 mm	5-0	\$113.99
Y494G	ETHICON	P-3	13 mm	4-0	\$381.62
Y494C	CP MEDICAL	NP-3	13 mm	4-0	\$111.02
Y495C	CP MEDICAL	NPS-2	19 mm	5-0	\$103.10
Y496G	ETHICON	PS-2	19 mm	4-0	\$356.82
Y496C	CP MEDICAL	NPS-2	19 mm	4-0	\$109.77
Y497G	ETHICON	PS-2	19 mm	3-0	\$353.27
Y497C	CP MEDICAL	NPS-2	19 mm	3-0	\$109.15

18" Violet

Y463C	CP MEDICAL	NP-3	18 mm	5-0	\$111.09
Y463G	ETHICON	P-3	18 mm	5-0	\$387.53
Y464C	CP MEDICAL	NP-3	18 mm	4-0	\$103.10
Y464G	ETHICON	P-3	18 mm	4-0	\$381.62

Vicryl - Polyglactin 910

18" Undyed Braided

⊖ J488G	ETHICON	P-1	11 mm	7-0	\$285.92
490C	CP MEDICAL	NP-1	11 mm	5-0	\$116.67
J492G	ETHICON	P-3	13 mm	6-0	\$273.87
492C	CP MEDICAL	NP-3	13 mm	6-0	\$104.72
J493G	ETHICON	P-3	13 mm	5-0	\$271.93
493C	CP MEDICAL	NP-3	13 mm	5-0	\$114.68
J494G	ETHICON	P-3	13 mm	4-0	\$269.35
494C	CP MEDICAL	NP-3	13 mm	4-0	\$114.68
J495G	ETHICON	PS-2	19 mm	5-0	\$251.65
495C	CP MEDICAL	NPS-2	19 mm	5-0	\$114.95
J496G	ETHICON	PS-2	19 mm	4-0	\$245.48
✓ 496C	CP MEDICAL	NPS-2	19 mm	4-0	\$117.01
J497G	ETHICON	PS-2	19 mm	3-0	\$243.59
497C	CP MEDICAL	NPS-2	19 mm	3-0	\$112.13
⊖ J499G	ETHICON	PS-3	16 mm	6-0	\$255.21
J507G	ETHICON	PS-4	16 mm	4-0	\$249.30

18" Violet Braided

J464G	ETHICON	P-3	13 mm	4-0	\$272.92
464C	CP MEDICAL	P-3	13 mm	4-0	\$113.32
⊖ J392H*	ETHICON	FS-2	19 mm	4-0	\$404.07
392C	CP MEDICAL	FS-2	19 mm	4-0	\$70.76
393C	ETHICON	FS-2	19 mm	3-0	\$70.76

27" Undyed Braided

⊖ VCP428H*	ETHICON	PS-2	19 mm	2-0	\$818.05
------------	---------	------	-------	-----	----------

30" Undyed Braided

422C	CP MEDICAL	FS-2	19 mm	4-0	\$75.68
423C	CP MEDICAL	FS-2	19 mm	3-0	\$68.84
B945C	CP MEDICAL	CP-1	37 mm	2-0	\$54.60

FA Sutures are fast absorbing.

*Sutures with item numbers ending in H are packaged in boxes of 36. All other sutures are sold in boxes of 12.

ABSORBABLE Polydioxanone 18" Clear Mono

Item	Make	Needle Style	Needle Size	Suture Size	Price
Z494G	ETHICON	P-3	13 mm	4-0	\$347.38
Z496G	ETHICON	PS-2	19 mm	4-0	\$369.82

30" Violet

M397C	CP MEDICAL	FS-2	19 mm	4-0	\$72.49
M452C	CP MEDICAL	FS-1	24 mm	3-0	\$76.21

Raising the Standard



Featuring

Swann-Morton®

✓ Consistency

Long lasting edge with tactile feel for better control

✓ ISO Certified

Blades and handles designed to fit for better control

✓ Reliability

Zero blade defects, no coating needed

✓ All-Inclusive

Product portfolio of blades, handles, disposable scalpels, miniature blades, and safety items available



Scan to view the Swann-Morton Collection

100 Blades Per Box*

*standard blade sizes





AccuThrive DermaBlade®

- MicroCoat® Coated Stainless Steel blade reduces skin and tissue trauma
- Flexible design for maximum control
- Single-use
- Sterile
- CE marked

Item	Description	Price
<input checked="" type="checkbox"/> DB	50/box	\$160.95

Dermaplaning Blade

Ideal for manual exfoliation of the epidermis

- Butter blade style
- Stainless steel blade
- Swann-Morton blades features an open edge grinding technique that retains sharpness
- Southmedic blades are polymer coated

Item	Description	Price
01SM10R	Swann-Morton, #10R, 100/box	\$59.04
<input checked="" type="checkbox"/> 730410R	Southmedic, #10R, 50/box	\$92.05

Southmedic Plus Blades

- Polymer-coated stainless steel
- Sterile and disposable
- CE marked

Item	Description	Price
730410	#10, 50/box	\$42.22
730411	#11, 50/box	\$42.22
730415	#15, 50/box	\$42.22
734415	#15C, 50/box	\$117.99

Swann-Morton® Blades

- Open edge grinding technique that retains sharpness
- Stainless steel
- Sterile and disposable
- CE marked

Item	Description	Price
#10/CS	#10, 100/box	\$50.00
#11/CS	#11, 100/box	\$50.00
#15/CS	#15, 100/box	\$50.00
#15C/CS	#15C, 100/box	\$97.00

Swann-Morton® Carbon Steel Blades

- Open edge grinding technique that retains sharpness
- Carbon Steel
- Sterile
- CE Marked

Item	Description	Price
00SM11	#11, 100/box	\$50.00
00SM15	#15, 100/box	\$50.00



Enhance Safety and Accuracy with Polymer Coated Blades

Southmedic's Precision Plus Safety Scalpels and 10R Dermaplaning Blades deliver exceptional precision and safety. The PPSS provides smooth, controlled cuts with minimal tissue trauma, while the 10R blades ensure precise, non-invasive exfoliation. Both feature sterile stainless steel blades with individual packaging.



Precision Plus Safety Scalpels


- Polymer-Coated Edge for smoother cuts
- Enhanced Safety with retractable, self-locking guard

10R Dermaplaning Blades

- Visibly Sharp Blade for smooth dermaplaning
- Easy to Use with standard reusable handle



Whether performing delicate dermaplaning or making precise surgical incisions, trust Southmedic for the precision you can rely on.



AccuThrive Razor Blades

Item	Description	Price
ASR/1	Double Edge, 100/pack	\$76.53
ASR/2	Double Edge, 250/pack	\$160.02
ASR/6	Single Edge, 250/pack	\$229.60



Surgery – Scalpels & Blades



Lumohs LED Lighted Scalpel

Trust Lumohs for better visualization when cutting around critical nerves, vessels and target tissue. Lumohs LED lighted scalpel makes surgery easier, more efficient and safer with fewer dark spots and shadows from obstructions caused by conventional lighting. Helps reduce surgical strain and errors due to poor visualization.

- Enhanced ergonomics for elliptical, round and Mohs excisions
- Improves efficiency and beveled margins
- Portable for bedside consults, shave biopsies and routine surgical procedures
- Reusable/autoclavable

Item	Description	Price
NLS-1612-009	Lumohs LED Lighted Scalpel Handle w/Battery	\$195.00
NLS-1614-001	Replacement batteries, 8/pk	\$51.75



Disposable Blade Remover

- Holds 100+ blades
- Blade drops automatically into reservoir

Item	Description	Price
BR/D	Plastic, with Reservoir	\$2.58



Quicksmart® Blade Remover

- Single-hand blade remover and wall mount bracket
- Holds 100 blades

Item	Description	Price
48QS	3¾" W x 6" H	\$38.80
48QSB	Wall Mount	\$9.00



DeRoyal® SurgiMate® Skin Staplers

Regular staple dimensions after closure: 3.7mm H x 5.7mm W. Wide staple dimensions after closure: 4.7mm H x 6.5mm W.

- Clear viewing window for maximum visibility
- 35 staples per stapler
- "Pre-cocking" feature for precision placement
- 316L stainless steel staples

Item	Description	Price
25-3002-6	Regular, Sterile, 6/box	\$170.21
25-3001-6	Wide, Sterile, 6/box	\$191.33



Delasco Blade Remover

- Stainless steel
- Autoclavable

Item	Description	Price
BR-101	Each	\$25.88

SurgiSeal Stylus® Topical Skin Adhesive

- Substantially equivalent to Dermabond® at a significantly lower price
- Rapid set time creating a better patient experience
- Microbial barrier - only product on the market with antimicrobial properties
- Greater service coverage when compared to other skin adhesives
- A single, thick, continuous layer may be sufficient to hold approximated wound edges closed for 5 to 10 days



Item	Description	Price
SS-050S	0.50 ml Applicator, 12/box	\$255.74



DERMarker™

- One time use, EZ Removable Ink®
- White ink visible on all skin tones
- Red ink preferred during laser procedures

Item	Description	Price
✓ 1444/G-50	Green, 50/pack	\$33.48
1444/R-50	Red, 50/pack	\$33.49
1444/W-50	White, 50/pack	\$37.43

Delasco Marking Ink

- Sometimes called Bonnie's Blue Ink
- Used to make persistent marks for outlining surgical procedures, radiation therapy areas, or patch tests
- Contains Gentian Violet (crystal violet), brilliant green, alcohol and water



Item	Description	Price
INK/1	1 oz	\$15.34
INK/2	2 oz	\$21.65 ◊

Universal Facial Wrap

- Includes a 3 1/2" soft segmented elastic wrap with contact closure which provides a soft, comfortable, contoured fit on the chin, cheeks or ears
- To prevent the band from sliding off, it also comes with a 1 1/2" soft segmented elastic stabilizing band with contact closure.
- Elastic face wrap 1 1/2" stabilizing band with contact closure
- Ideal for otoplasty and liposuction of the chin and cheeks
- Universal size
- Latex-free

Item	Description	Price
2020	3 1/2" wide	\$18.54



Facial Garment

- Compression for chin and neck regions
- No internal seams
- Latex-free

For sizing, measure head circumference above ears. Sizes available for 8006: (S)mall 19.5"-20.5", (M)edium 20.5"-21.5", (L)arge 21.5"-22.5", (XL)arge 22.5"-23.5" and (XXL)arge 23.5"-24.5"

Item	Description	Price
8006M	Single closure, medium	\$16.22
8006L	Single closure, large	\$18.78
8006XL	Single closure, x-large	\$18.78
8006XXL	Single closure, xx-large	\$16.22
2002	One Size, 5 Closures	\$30.50



Dermoscopy – DermLite® Dermatoscopes & Lumio® Viewers

DermLite® Dermatoscopes

Every DermLite dermatoscope comes with a characteristic color spectrum designed to highlight features in skin lesions that are critically important in early detection. Over the past 20 years, DermLite has perfected that blend of colors, another reason why it represents best-in-class in dermoscopy.

- Mobile compatible and rechargeable battery except DL100 model
- Used in viewing pigmented, vascular, or inflammatory lesions, varicose veins, skin irritants, and fungal infections



D-LITE/200HR

Item	Model	Price	Optics	Immersion Fluid	C/W Camera	10X Mag	10mm Reticle	Polarized	Includes
D-LITE	DL100	\$395.00	15 mm	-	-	•	•	Cross	B,P,I,S
⊖ D-LITE/1	DL1	\$525.00	15 mm	Optional	•	•	•	Cross & Non	C
D-LITE/200HR	HR 200	\$1,050.00	25 mm	-	-	•	•	Cross	B,C,P,I,S
✓ D-LITE/200HYBRID	HYBRID 200	\$1,050.00	25 mm	Optional	•	•	•	Cross & Non	B,C,P,I,S
✓ D-LITE/200HYBRID/G	(Gold)	\$1,050.00	25 mm	Optional	•	•	•	Cross & Non	B,C,P,I,S
D-LITE/3N	DL3N	\$1,395.00	25 mm	Optional	•	•	•	Cross & Non	B,C,P
D-LITE/4W	DL4W	\$1,595.00	30 mm	Optional	•	•	•	Cross & Non	B,C,P,I,S
✓ D-LITE/4	DL4	\$1,595.00	30 mm	Optional	•	•	•	Cross & Non	B,C,P,I,S
✓ D-LITE/5	DL5	\$1,795.00	32 mm	-	•	6X/1.4X	•	Cross & Non	B, C, P, I

*Includes: Battery (B), Charger (C), Belt Pouch (P), Ice Cap Covers (I), Silicone Sleeve (S)

Let our experts help you select the right dermatoscope. Contact our sales team at 800.831.6273.



Lumio® 2 Hand Held Cross Polarized Viewer with UV

By combining brightest-in-class polarized white light, three wavelengths in the UV spectrum, a special Woods mode, and special optical filters in a beautifully slim, quickly rechargeable device, Lumio 2 is the ultimate tool for visualizing various skin conditions.

Includes one Lumio 2, OptiClip 405nm Long-Pass, OptiClip 2.5x, USB-C to USB cable, and one protective neoprene pouch.

- 100mm 2x lens
- White light at 365nm, 385nm, 405nm and Woods
- OptiClip filter for enhanced fluorescence imaging
- OptiClip factor 2.5 magnification boost



Item	Description	Price
✓ LUMIO2	2x Magnification, UV light only	\$995.00

Dermlite® Dermatoscopes & Lumio® Viewers – Dermoscopy

Lumio®

- Includes device, four AA batteries and one protective neoprene pouch.
- 40 bright-white LEDs for LUMIO
 - 40 UV LEDs for LUMIO/UV
 - Ideal for general skin exams
 - 75 mm lens
 - 5-year warranty



Lumio® S

- Includes USB charging cable.
- 40 mm lens
 - Cross-polarized and non-polarized
 - Large field of view
 - 5-year warranty



Item	Description	Price	Item	Description	Price
<input checked="" type="checkbox"/> LUMIO	2X Magnification, white light	\$595.00	LUMIO/S	4x Magnification	\$625.00
LUMIO/UV	2x Magnification, UV light	\$595.00	<input type="checkbox"/> LUMIO/S/SK	Silicone Sleeve	\$25.00

Dermlite DL5 with UV

Become ultra visionary.



©2024 All rights reserved. Dermlite is a registered trademark | See dermlite.com/patents for more details. Tinea image ©Kasia Kordecka





Duo-Check™ Pouch

- Built-in internal/external process indicators that change color when exposed to a specific temperature
- For use with steam or EO sterilization processes only
- Self-sealing

Item	Description	Price
CROSCXS	3½" x 5¼", 200/box	\$11.84
CROSCQ	3¼" x 9", 200/box	\$14.58
CROSCM	5¼" x 10", 200/box	\$27.53
CROSCCL	7½" x 13", 200/box	\$51.07
CROSCCL816	8" x 16", 200/box	\$71.81
Sure-Check Pouch (For steam and EtO only)		
CROSCS2	3½" x 9", 200/box	\$18.71



Dynarex Sterilization Pouch

- Built-in internal/external process indicators that change color when exposed to a specific temperature
- For use with chemical vapor, steam or EO sterilization processes only
- Self-sealing

Item	Description	Price
DYN4462	3½" x 9", 200/box	\$8.00
DYN4464	5¼" x 10", 200/box	\$13.44
DYN4465	7½" x 13", 200/box	\$24.00



EmPower™ Dual Enzymatic Detergent

- Dual-enzymatic detergent cleaner
- Gentle on instruments
- For use with automated or manual instrument processing

Item	Description	Price
MET10-4100	Gallon	\$55.43



Halyard CSR Wrap

- Tear, flame, lint and abrasion-resistant
- Enhanced microbial barrier protection
- Pre-vacuum steam sterilant penetration
- Low linting

Item	Description	Price
KIM37050CS	20" x 20", 100/pk, 10pk/case	\$249.35
KIM37051CS	24" x 24", 100/pk, 5pk/case	\$155.87

Curettes, Extractors & Elevators – Instruments



Dermal Curettes

Fox Curette, Oval 5½"

FC-2	2 x 3 mm	\$106.61
FC-3	3 x 4 mm	\$106.61
FC-4	4 x 6 mm	\$106.61
FC-5	5 x 7 mm	\$106.61
FC-6	6 x 9 mm	\$102.47

Fox Curette, Round 5½"

✓ FCR-2	2 mm	\$95.22
FCR-3	3 mm	\$95.22
✓ FCR-4	4 mm	\$95.22
FCR-5	5 mm	\$95.22
FCR-6	6 mm	\$95.22

Bone Curettes

Piffard Oval Curette 5½" (Miltex)

† 33-16-1	Narrow Handle, 3 x 3.5 mm	\$142.39
-----------	------------------------------	----------

Martini 5½" (Miltex)

21-321	Double End, 4 mm and 5 mm Round Cups	\$120.20
--------	---	----------

Williger 5½" (Miltex)

21-322	Double End, 3 mm and 4 mm Oval Cups	\$120.20
--------	--	----------

Comedone Extractors

Schamberg 3¾"

CE-11	Regular, 2 Pinched Ends	\$26.00
CE-14	Extra Fine, 2 Pinched Ends	\$27.64

CE-12	Regular, 1 Pinched & 1 Squared End	\$26.00
CE-13	Extra Fine, 1 Pinched & 1 Squared End	\$26.00
✓ CE-21	Regular, Square Ends	\$27.64
CE-23	Extra Fine, Square Ends	\$26.00

Comedone Extractor 4¼"

✓ CE-31	Flat End & Wire End Expressors	\$60.49
---------	-----------------------------------	---------

Saalfield

CE-54	4¼", w/Acne Lancet	\$26.00
CE-56	5¾", w/Acne Lancet	\$31.75

Delasco 6"

CE-72	Fenestrated Cup and One Delicate Curved Lancet w/Cap	\$52.79
-------	--	---------

Elevators

Freer Septum Elevator 7½"


E-122	Double End, Sharp and Blunt Blades, 5 mm Wide	\$59.15
-------	--	---------

Miltex Locke Elevator

40-94	4½" with 3 x 15 mm Blade	\$55.55
40-96	5" with 5 x 20 mm Blade	\$55.53

†May take 2-4 weeks to ship.

Instruments - Forceps



Tissue Forceps

Adson 4 3/4" Forceps

<input checked="" type="checkbox"/>	FA-1	Serrated, 1.5 mm Tips	\$20.70
	FA-5	Serrated, Fine, 1.0 mm Tips	\$34.16
	FA-1F	Serrated, Fenestrated Handle, 1.5 mm Tips	\$91.08
<input checked="" type="checkbox"/>	FA-2	1 x 2 Teeth, 1.5 mm Tips	\$22.77
	FA-6	1 x 2 Teeth, Fine, 1.0 mm Tips	\$35.19
	FA-2S	1 x 2 Teeth & Serrated, 1.5 mm Tips	\$57.16
	FA-2F	1 x 2 Teeth, Fenestrated Handle, 1.5 mm Tips	\$94.19
	FA-3	2 x 3 Teeth, 2.0 mm Tips	\$48.65
	FA-SPF	1 x 2 Teeth, w/Suture Tying Platform, Fenestrated Handle, 1.5 mm Tips	\$95.22
	FA-SP	1 x 2 Teeth, w/Suture Tying Platform, 1.5 mm Tips	\$84.87
	FA-9SP	1 x 2 Teeth, w/Suture Tying Platform, TC, 1.5 mm Tips	\$134.55

Brown-Adson 4 3/4" Forceps

	FA-4	7 x 7 Side Grasping Teeth, 2.0 mm Tips	\$57.57
--	-------------	--	---------

Bishop-Harmon 3 3/4" Forceps

	BH-1	Serrated, Fenestrated, 0.8 mm Tips	\$151.11
<input checked="" type="checkbox"/>	BH-2	1 x 2 Teeth, Fenestrated, 0.8 mm Tips	\$180.09
	BH-2SP	1 x 2 Teeth, Fenestrated, w/Suture Tying Platform, 0.8 mm Tips	\$194.58

Graefe Iris 2 3/4"


	GI-3	Serrated, Straight, Fine, 1.0 mm Tips	\$97.29
<input checked="" type="checkbox"/>	GI-7	1 x 2 Teeth, Straight, Fine, 1.0 mm Tips	\$105.57

Semken Iris 4", with Aligning Pin

	G-10	1 x 2 Teeth, Straight, Standard, 1.0 mm Tips	\$61.07
--	-------------	--	---------

*Denotes item pictured.

Forceps & Needle Holders – Instruments



Dissecting 5"

G-23	Extra Narrow, Serrated, 1.5 mm Tips	\$19.67
G-24	1 x 2 Teeth, Extra Narrow, 1.5 mm Tips	\$22.77

Splinter, Jewelers & Celia Forceps

Feilchenfeld Splinter		
SF-11	3" Serrated	\$21.74

Jewelers/Micrograft Forceps

J-3	4 3/4" Straight, Narrow, Fine, Style 3	\$25.88
J-5	4 3/8" Straight, Super Fine, Style 5	\$25.88
J-7	4 3/8" Straight, Super Fine, Style 5	\$37.26

Cilia Forceps

CF-1	3 1/2" Littauer, 4 mm Wide Jaws w/Fine Horizontal Serrations	\$19.67
-------------	--	----------------

CF-3	3 1/2" Bergh, 5 mm Wide Jaws w/Fine Horizontal Serrations	\$19.67
<input checked="" type="checkbox"/> CF-11	3 3/4" 3 mm Diagonal Smooth Jaw	\$30.17

Carmalt Splinter 4 1/4"

<input checked="" type="checkbox"/> SF-15	Straight, Serrated Tip w/Aligning Pin	\$43.77
--	--	----------------

Plain Splinter

SF-17	3 1/2" Serrated w/Aligning Pin	\$20.70
--------------	--------------------------------	----------------

Needle Holders

Webster 5"

<input checked="" type="checkbox"/> W-13	Smooth	\$43.47
<input checked="" type="checkbox"/> W-13S	Serrated Jaws	\$66.24
<input checked="" type="checkbox"/> W-13C	Smooth, Tungsten Carbide	\$115.92
W-13CS	Serrated, Tungsten Carbide	\$137.66

Webster 4 3/4"

W-12	Smooth, Delicate	\$94.27
W-12S	Extra Fine Serrations, Delicate	\$97.00

Instruments – Needle Holders, Clamps & Hemostats



Halsey 5"		
<input checked="" type="checkbox"/>	WH-12 Smooth Jaws	\$57.85
	WHS-12 Fine Serrations	\$49.06
	WH-93 Fine Serrations, Tungsten Carbide	\$130.44
Baumgartner 5½"		
	B-13CS Serrated Jaws, Tungsten Carbide	\$108.68
	B-13S Serrated Jaws	\$54.86
Derf 4¾"		
<input checked="" type="checkbox"/>	D-12C Smooth Jaws, Tungsten Carbide	\$151.11
Needle Holder / Suture Scissors Combination		
Olsen-Hegar 4¾"		
	OHS-12C Serrated, Fine Jaws, Tungsten Carbide	\$185.27

Towel Clamps

Backhaus

<input checked="" type="checkbox"/>	TC-09 3½"	\$33.12
	TC-11 4"	\$34.16

Hemostats

Pedifine Hartman Mosquito 3½"

	H-112 Curved, Extra Fine	\$36.23
--	---------------------------------	---------

Hartman Mosquito 3¾"

	H-116 1 x 2 Teeth, Curved	\$65.21
--	----------------------------------	---------

Halsted 3¾"

	H-152 Curved	\$35.19
--	---------------------	---------

Micro-Mosquito 4¾"







	H-122 Curved, Fine	\$40.37
--	---------------------------	---------

	H-123 Straight, Fine	\$39.33
--	-----------------------------	---------

Halsted Mosquito 5" Hemostats

<input checked="" type="checkbox"/>	H-140 Curved, Standard	\$50.64
	H-141 Straight, Standard	\$48.32

Hemostats & Nail Nippers – Instruments



Kelly 5½"

H-182	Curved	\$34.16
H-183	Straight	\$33.12

Providence Hospital 5½" (Miltex)

7-30	Straight, Serrated Jaws	\$97.45
------	-------------------------	---------

Nail Nippers

Nail Nipper 5"

<input checked="" type="checkbox"/>	N-26	Straight Jaws, Double Spring, Stainless	\$127.19
-------------------------------------	------	---	----------

Nail Nipper 6"

N-27	Straight Jaws, Double Spring, Stainless Steel	\$149.04
------	---	----------

Dual Action 6"

<input checked="" type="checkbox"/>	N-7	Straight Jaws, Sharp Stainless Steel	\$512.33
-------------------------------------	-----	--------------------------------------	----------

Barrel Spring 5½"

N-3	Concave Jaws, Stainless Steel	\$118.90
-----	-------------------------------	----------

Economy Nail Nippers 5½"

<input checked="" type="checkbox"/>	DNN-1	Concave, Stainless Steel	\$15.69
<input checked="" type="checkbox"/>	DNN-2	Straight, Double Spring, Stainless Steel	\$15.69

Cuticle/Nail Nippers

N-6	Straight Jaws, Double Spring, Fine, Stainless Steel	\$112.74
-----	---	----------

V940200	Concave, Single Spring, Stainless (Miltex Vantage)	\$50.74
---------	--	---------

Instruments – Disposable




DK-1PK



DK-2PK



DK-3PK




SS62



SS63



SS65



SS98



SS66



SS69



SS67



SS68

Disposable Instrument Kits

Biopsy Kit

Includes Iris scissors (DS-8), needle holder (DNH-1), Adson forceps (DF-4) and 5¼" X 10" sealable transparent autoclavable pouch.

<input checked="" type="checkbox"/>	DK-1PK	Biopsy Kit, 12/set	\$81.79
	DK-2PK	Biopsy Kit with Punch, 12/set	\$117.64

Excision Kit

Includes Metzenbaum scissors (DS-4), Iris scissors (DS-8), needle holder (DNH-1), Adson forceps (DF-4), hemostat (DH-2), #3 scalpel handle and 5¼" x 10" sealable transparent autoclavable pouch.

DK-3PK	Excision Kit, 6/set	\$102.84
---------------	---------------------	-----------------

Sterile Skin Biopsy Kit

Includes: 30G x ½" needle, 18G x 1" needle, three PVP swabsticks, three alcohol swabsticks, CSR wrap, tray, Adson 1 x 2 teeth forceps, #11 mini scalpel, 3 mm biopsy punch, six 2" x 2" gauze, 1 cc syringe and fenestrated drape.

DK-S1	Skin Biopsy Kit	\$18.84
--------------	-----------------	----------------

Sterile Suture Removal Kit

Includes 4½" Littauer scissors, 4" thumb forceps and cotton filled sponge in sterile tray package. 10 disposable sterile kits per pack.

101L	10 Kits	\$32.02
-------------	---------	----------------

Sterile Instruments

Medium Grade

Disposable. Packages of 10 sterile individually wrapped. Made in Pakistan.

SS62	Hemostat, Mosquito, 5" Straight, 10/pack	\$30.75
SS63	Hemostat, Mosquito, 5" Curved, 10/pack	\$30.96
SS65	Needleholder, Webster, 5" Smooth, 10/pack	\$31.56
SS98	Needleholder, Halsey, 5" Serrated, 10/pack	\$32.98
SS66	Scissor, Iris, 4½" Straight, 10/pack	\$30.36
SS69	Scissor, Iris, 4½" Curved, 10/pack	\$26.68
SS67	Forceps, Adson 4¾", 1 X 2 Teeth, 10/pack	\$23.45
SS68	Forceps, Adson 4¾", Serrated, 10/pack	\$22.76

Disposable & Dextra by Delasco – Instruments

Scissors

DS-7	Iris, 4½" Straight	\$2.25
✓ DS-8	Iris, 4½" Curved	\$2.25

Forceps

✓ DF-1D	4½" Tissue, 1 x 2 Teeth	\$2.25
✓ DF-3	4¾" Adson, Serrated	\$2.25
✓ DF-4	4¾" Adson, 1 x 2 Teeth	\$2.25

Scalpel Handles

DSH	#3 Handle, Stainless, 5"	\$2.25
------------	--------------------------	--------

Needle Holders

✓ DNH-1	5" Webster, Smooth Jaws	\$4.48
DNH-4	4¾" Derf	\$4.48

Comedone Extractors

DCE-1	3¾" Regular, Square, Schamberg	\$3.36
DCE-2	4" Extra-Fine, One Pinched, Schamberg	\$4.48



DEXTRA BY DELASCO INSTRUMENTS

Reusable instruments, but economical enough to be disposed when working on high-risk patients. Dextra by Delasco offers superior quality at an unbelievable price.

Dermal Curettes

FCR-3M	Fox 5½", 3 mm	\$14.49
FCR-4M	Fox 5½", 4 mm	\$14.49

Comedone Extractors

CE-21M	Schamberg 3¾", Regular, Square Ends	\$12.42
CE-31M	4¼", Flat End & Wire End Expressors	\$12.42
CE-72M	Hallack 6", Fenestrated Cup	\$12.42

Scissors, Fine

S-117M	Iris 4½" Straight, S/S	\$9.32
S-118M	Iris 4½" Curved, S/S	\$9.32
S-928M	Iris 4½" Curved, S/S, Tungsten Carbide	\$18.63
S-702M	Gradle 3¾" Slightly Curved, S/S, Fine	\$9.32

Scissors, Operating

S-562M	Metzenbaum 5" Curved, B/B, Fine	\$12.42
S-563M	Metzenbaum 5" Straight, B/B, Fine	\$12.42

Bandage Removal

S-007M	Lister 5¾"	\$14.49
---------------	------------	---------

Suture Removal

S-011M	Spencer 3½" Satin Finish, Fine	\$9.32
---------------	--------------------------------	--------

Skin Hooks

SH-212M	Frazier 5" Sharp	\$9.32
----------------	------------------	--------

Scalpel Handles

H-323DM	#3 - Fits Blades #10-#15C, 5" w/Metric Rule	\$8.28
----------------	---	--------

H-473M	Round - Fits Blades #10-#15C, 6" Straight, Knurled	\$12.42
---------------	--	---------

Forceps

FA-1M	Adson 4¾", Serrated Jaws	\$9.32
FA-2M	Adson 4¾", 1 x 2 Teeth	\$9.32
FA-SPM	Adson 4¾", 1 x 2 Teeth w/Tying Platform	\$12.42
FA-6M	Micro-Adson 4¾", 1 x 2 Teeth	\$12.42
BH-1M	Bishop-Harmon 3¾", Serrated, Fenestrated, 0.8 mm Tips	\$14.49
BH-2M	Bishop-Harmon 3¾", 1 x 2 Teeth, Fenestrated, 0.8 mm Tips	\$14.49

Needle Holders

WH-12M	Halsey 5" Smooth Jaws	\$14.49
W-13CM	Webster 5" Smooth Jaws, Tungsten Carbide	\$22.77
W-13M	Webster 5" Smooth Jaws	\$14.49

Towel Clamps

TC-09M	Backhaus 3½"	\$14.49
---------------	--------------	---------

Hemostats

H-140M	Halsted Mosquito 5" Standard, Curved	\$12.42
---------------	--------------------------------------	---------

Nail Nippers

N-7M	Nipper 6", Straight Jaws	\$138.69
N-32M	Anvil Style Splitter 5", Double Spring	\$45.54

Instruments – Operating Scissors & Suture



Scissors, Operating

Lexer-Baby 4"

S-110 Curved, B/B \$52.24

Plastic Surgery, Wagner 4¾"

S-811 Straight, S/S, 30 mm Blades \$42.44

S-812 Curved, S/S, 30 mm Blades \$43.47

S-912 Curved, S/S, 25 mm Blades, Tungsten Carbide \$174.30

Metzenbaum

S-562 5" Curved, Fine, B/B, 25 mm Blades \$44.51

S-563 5" Straight, Fine, B/B \$43.47

S-572 5¾" Curved, B/B, Standard, 35 mm Blades \$49.68

S-574 5¾" Curved, Fine S/S, 35 mm Blades \$65.21

Mini Metzenbaum 4½"

S-558 Curved, B/B, 25 mm Blades \$46.58

Kilner 4½"

S-938 Curved, Flat Tip, B/B, Serrated, Tungsten Carbide, Large Rings \$214.77

Ragnell

S-322 4¾" Slight Curve, B/B, Flat Serrated Tip \$78.60

Deaver 5½"

S-021 Curved, S/S \$52.79

S-031 Straight, S/B, 40 mm Blades \$41.40

Suture Applying

Oliver 4½"

S-101 Straight, S/S Tips (Delasco Exclusive) \$35.19

S-103 Curved S/S Tips (Delasco Exclusive) \$37.26

Suture Removal

Spencer

S-011 3½" Satin Finish, Fine \$40.37

S-010 4½" Fine \$77.13

Suture, Scalpel Handles & Skin Hooks - Instruments



Shortbent 3½"

S-013 Curved Upwards, Fine \$77.01

Littauer

S-014 5" Light Pattern \$38.30

Buck 5¾"

S-016 Straight \$69.38

O'Brien 3¾"

S-019 Angled Sharp Points \$88.52

O'Brien/Spencer 4½"

S-020 Angled, with Spencer Style Blade \$85.06

Scalpel Handles

#3 Handles - Fits Blades #10-#15C

✓ H-323D 5" with Metric Rule \$21.26

H-333 5" Safety Scalpel Handle \$34.81

Round Handles - Fits Blades #10-#15C

H-473 6" Straight, 8.25 mm, Round Handle \$50.06

BR-101 Removes Blades from #3 and #7 Handles \$25.88

Skin Hooks

O'Connor 4¾"

SH-105 3 mm Diameter, Sharp \$94.19

Frazier 5"

✓ SH-212 3.5 mm Diameter, Sharp \$91.08


SH-213 3.5 mm Diameter, Blunt \$91.08

Guthrie

SH-162 4¾" 1 mm Diameter, 2 Sharp Prongs, 1.5 mm Wide \$130.41

SH-172 5" 1.5 mm Diameter, 2 Sharp Prongs, 2 mm Wide \$130.41

Instruments – Scissors, Fine



Scissors, Fine Gradle

- ✓ **S-702** 3¾" Curved, S/S, 15 mm Fine Blades \$94.53

Iris 4½" Supercut

One blade honed to razor edge. Marked for supercut identification

- S-119** Straight, S/B, 20 mm Blades \$64.60
- S-120** Curved, S/B, 20 mm Blades \$95.97
- ✓ **S-137** Curved, S/S, 25 mm Blades, Serrated \$175.95
- S-137B** Curved, S/S, 25 mm Blades \$178.02
- S-139** Straight, S/S, 25 mm Blades, Serrated \$176.99
- S-139B** Straight, S/S, 25 mm Blades \$171.81

Iris

- ✓ **S-117** 4½" Straight, S/S, 25 mm Blades \$38.30
- ✓ **S-118** 4½" Curved, S/S, 25 mm Blades \$40.37
- S-318** 4½" Curved, S/S, 25 mm Blades, Serrated \$43.47
- S-126** 4½" Curved, S/B, 20 mm Blades \$62.64
- S-146** 4½" Curved, B/B, 20 mm Blades \$73.62

- S-927** 4½" Straight, S/S, 25 mm Blades, Tungsten Carbide \$120.06
- S-927SC** 4½" Straight, S/S, 25 mm Blades, Serrated, Supercut, Tungsten Carbide \$232.88
- ✓ **S-928** 4½" Curved, S/S, 25 mm Blades, Tungsten Carbide \$124.20
- S-928SC** 4½" Curved, S/S, 25 mm Blades, Serrated, Supercut, Tungsten Carbide \$239.09
- S-114** 4" Curved, S/S, 20 mm Blades \$49.68
- S-115** 4" Straight, S/S, 20 mm Blades \$47.61
- S-112** 3½" Curved, S/S, 20 mm Blades \$45.54
- S-113** 3½" Straight, S/S, 20 mm Blades \$43.47

Microscissors

- MS-15** 6" Castroviejo, Curved, S/S, 12 mm Blades, Serrated \$670.68

*Denotes item pictured.

Scissors, Fine & Bandage Removal – Instruments



LaGrange 4½"

S-711 Compound Curve, S/S, Serrated, For Undermining Dome Shapes or Clipping of Donor Plugs **\$75.56**

Stevens Tenotomy

S-721 4¼" Curved, S/S, 20 mm Blades **\$50.72**

S-730 4¼" Straight, B/B, 20 mm Blades **\$45.54**

Scissors, Bandage Removal

Lister

S-005 7¼" **\$52.77**

S-006 4¾" **\$35.19**

S-007 5¾" **\$39.33**

5-550 8" with Larger Finger Ring (Miltex) **\$162.23**

Nurse's Bandage Removal (Miltex)

Fluoride coated for easier clean-up.

5-700 5½" Pink Plastic Finger Ring & Safety Guard **\$52.92**

5-702 5½" Blue Plastic Finger Ring & Safety Guard **\$52.92**


Utility Scissors (Miltex Vantage)

V95-1026 7½" Orange Handle **\$10.12**

Scissors, Undermining (Miltex)

33-242 4½" Curved, Flat Pointed Tips **\$135.52**

For a full line of instruments go to delasco.com



TSK STERiGLIDE™

Aesthetic Cannulas

The STERiGLIDE™ outperforms any other cannula available and leads the market as the gold standard. Ultra-thin wall made of surgical stainless steel with a polycarbonate hub.

- Dome-shaped tip
- Most distal port location of all leading brands
- Proprietary surface treatment
- Hub clearly marked with port position
- Sterile, CE marked

Each cannula packaged with a sharp needle of optimal gauge for entry hole.

Item	Gauge & Size	Price
2238SG	22G x 1½", 20/box	\$184.53
2250SG	22G x 2", 20/box	\$184.53
2270SG	22G x 2¾", 20/box	\$184.53
✓ 2538SG	25G x 1½", 20/box	\$184.53
2550SG	25G x 2", 20/box	\$184.53
2725SG	27G x 1", 20/box	\$184.53
2738SG	27G x 1½", 20/box	\$184.53
2750SG	27G x 2", 20/box	\$184.53
3025SG	30G x 1", 20/box	\$184.53



TSK CSH Cannulas

- Ultra-thin wall
- Blunt round end
- Sterile, CE marked

Item	Gauge & Size	Price
2250TSK	22G x 2", 25/box	\$138.31
2538TSK	25G x 1½", 25/box	\$138.31
2550TSK	25G x 2", 25/box	\$138.31
2738TSK	27G x 1½", 25/box	\$138.31
2750TSK	27G x 2", 25/box	\$138.31




DermaSculpt Microcannulas

- Blunt round end, sterile

Includes sharp needles of a corresponding gauge.

Item	Gauge & Size	Price
1870DS	18G x 2¾", 20/box	\$117.00
2250DS	22G x 2", 20/box	\$117.00
2329DS	23G x 1⅛", 20/box	\$117.00
NEW 2538DS	25G x 1½", 20/box	\$117.00
✓ 2550DS	25G x 2", 20/box	\$117.00
2725DS	27G x 1", 20/box	\$117.00
2738DS	27G x 1½", 20/box	\$117.00
2750DS	27G x 2", 20/box	\$117.00
3025DS	30G x 1", 20/box	\$117.00



Air-Tite Aesthetics

THE FIRST 0.5 mL LUER LOCK LOW DEAD SPACE SYRINGE



SHOP NOW!



TSK Products

High Quality Specialty Hypodermics



STERIGLIDE™
PREMIUM MICROCANNULAS



STERIJECT™
HYPODERMIC NEEDLES

Luer Lock Syringe Connector


- Female to Female
- Disposable polypropylene

Item	Description	Price
DSA/S	Sterile, Green, 10/pack	\$23.53
11397PK	Non-sterile, White, 20/pack	\$16.81



Exel Safety-Lok™ Needle

Item	Gauge & Size	Price
EXE27437	30G x 1/2", 100/box	\$23.62




BEC305757

EXE27437

BD Eclipse Safety Needle

Item	Gauge & Size	Price
BEC305757	30G x 1/2", 100/box	\$43.57


Sharps – Needles



Exel Hypodermic Needle

- Regular bevel tip

Item	Gauge & Size	Price
EXE26439	30G x 1", 100/box	\$9.48
2834	28G x 3/4", 100/box	\$20.81
EXE26400	27G x 1/2", 100/box	\$8.77
EXE26427	27G x 1 1/4", 100/box	\$8.77
EXE26403	25G x 5/8", 100/box	\$8.77
EXE26405	25G x 1", 100/box	\$8.77
EXE26406	25G x 1 1/2", 100/box	\$8.77
EXE26412	22G x 1 1/2", 100/box	\$8.77
EXE26414	21G x 1", 100/box	\$8.77
EXE26416	21G x 1 1/2", 100/box	\$8.77
EXE26418	20G x 1 1/2", 100/box	\$8.77
EXE26419	18G x 1", 100/box	\$8.77



BD PrecisionGlide

- Regular bevel tip

Item	Gauge & Size	Price
3012/BD	30G x 1/2", 100/box	\$51.18
301/BD	30G x 1", 100/box	\$51.18
BEC305109	27G x 1/2", 100/box	\$18.06
BEC305122	25G x 5/8", 100/box	\$18.06
BEC305127	25G x 1 1/2", 100/box	\$18.06
305165	21G x 1", 100/box	\$18.34
305167	21G x 1 1/2", 100/box	\$18.07



TSK STERIJECT™

- Tri-beveled tip
- Siliconized shaft
- CE Marked

Item	Gauge & Size	Price
✓ 3312S	33G x 1/2", 100/box	\$66.02
✓ 3212S	32G x 1/2", 100/box	\$47.35
3112S	31G x 1/2", 100/box	\$42.44
3012/UTW	30G x 1/2", 100/box	\$55.18



Economy Hypodermic Needle

Item	Gauge & Size	Price
3012	30G x 1/2", Nipro, 100/box	\$11.37
3034	30G x 3/4", Covidien, 100/box	\$16.87



Vacutainer Safety-Lok™ Blood Collection

- Safety-engineered winged set for infusion and blood collection Luer Lock Tip
- Shield allows for one-handed activation
- Minimizes the risk of needlestick injuries
- Clear visibility of blood flashback

Item	Gauge & Size	Price
BEC367292	23G x 3/4", 7" tubing, 50/box	\$148.22
BEC367281	21G x 3/4", 12" Tubing, 50/pack	\$142.53
BEC364815	Vacutainer Disposable Holder, 250/pack	\$31.03



Sharps Containers

(From left to right)

Item	Description	Price
✓ DYN4624	Dynarex Mailbox, 5 Quart: 11.8" H x 10.7" W x 4.6" D	\$6.83
BEC305053	BD, 3 gallon: 15½"H x 13½"W x 6"D	\$18.54
MDG8708	Medegen Mailbox, 5.4 Quart: 10.5" x 3.5" x 10.5"	\$6.72
6087	Covidien, 2 gallon: 11"H x 10½"W x 6¾"D	\$9.25
KEN8507SA	Cardinal Mailbox, 5 Quart: 12½" H x 10¾" W x 5½" D	\$8.59

NeedleSafe II™

Single-handed needle uncapper, recapper and syringe holder. Hazard orange color. Base mounts horizontally or vertically to workplace with strong adhesive strip.

- NeedleSafe assembly slides off base and can be steam autoclaved, gas sterilized or disinfected by wiping with or soaking in alcohol.
- Universal gripper core inside assembly securely grips needle caps for quick, single-handed uncapping, recapping and safe disposal.
- Complies fully with OSHA regulations.
- NeedleSafe should be replaced every 4-6 months or immediately if the adhesive base becomes loose.

Sold individually.

Item	Price
NS2	\$25.88



0.5mL Luer Lock Syringe

- Easy to read 0.01mL graduations for precise dosing
- Low Dead Space plunger tip for minimal loss

Item	Gauge & Size	Price
MINI05	MinimLL 0.5mL Luer Lock Syringe	\$39.95

NEW



Sharps – Syringes



Ulticare™ Syringe with Fixed 31G Needle

Item	Gauge & Size	Price
<input checked="" type="checkbox"/> 3103/BT	0.3 ml, 31G x 8 mm, 100/box	\$23.30
31036/BT	0.3 ml, 31G x 6 mm, 100/box	\$23.30
3105/BT	0.5 ml, 31G x 8 mm, 100/box	\$23.30
31056/BT	0.5 ml, 31G x 6 mm, 100/box	\$23.30
311/BT	1 ml, 31G x 8 mm, 100/box	\$23.30

BD Syringe with Fixed 31G Needle

Item	Gauge & Size	Price
<input checked="" type="checkbox"/> 31516B	0.3 ml, 31G x 8 mm, 100/box	\$86.92
31036/BD	0.3 ml, 31G x 6 mm, 100/box	\$80.48
31516C	0.5 ml, 31G x 8 mm, 100/box	\$80.48
31056/BD	0.5 ml, 31G x 6 mm, 100/box	\$80.48
31516	1 ml, 31G x 8 mm, 100/box	\$82.65
31106/BD	1 ml, 31G x 6 mm, 100/box	\$75.15



BD Syringe/Needle Combo

- Luer Lock Tip
- 1 ml: 1/100 ml Graduations

Item	Gauge & Size	Price
BEC309626	1 ml, 25G x 5/8", 100/box	\$34.53
BEC309623	1 ml, 27G x 1/2", 100/box	\$38.34

1 ml Syringe

- Latex-Free
- 1/100 ml graduations

Item	Description	Price
92000/BD	BD Luer-Lok™, 100/box	\$90.92
92000	Exel Luer Lock, 100/box	\$23.54
92000/L	Exel Luer Lock, Low Dead Space, 100/box	\$28.87
EXE26048	Exel Luer Slip, Low Dead Space, 100/box	\$14.43
<input checked="" type="checkbox"/> BS/1	Norm-Ject Luer Slip, Low Dead Space, 100/box	\$19.95

3 ml Syringe

- Luer lock
- 1/10 ml graduations

BEC309657	BD, 200/box	\$22.71
-----------	-------------	---------

5 ml Syringe

- Luer lock
- 1/5 ml graduations

BEC309646	BD, 125/box	\$32.62
-----------	-------------	---------

10 ml Syringe

- Luer lock
- 1/5 ml graduations
- 1/20 ml graduations

LF/10	Nipro, 100/box	\$22.17
BEC309695	Control Tip, 25/box	\$67.04



Bovie® Smoke Shark™ III

Uniquely designed for the physician's office. The Bovie® Smoke Shark™ III is specifically intended for shorter (<15 minutes) procedures. Compact size, ultra-quiet (<55dB) and powerful vacuum for efficient smoke evacuation. Multi-procedure 4-Stage ULPA filter to remove and trap 0.1-0.2µm particles and micro-organisms at 99.999% efficiency. This unit also comes with a remote activation switch that starts the vacuum when the electrosurgery unit is activated and automatically stops when it's deactivated to maximize filter life.



- Includes: (1) 6-hour filter, SmartShark remote activation switch, (1) SharkSkin pencil adapter w/10' tubing
- 3/8" port filtration tubing connection
 - 9.4" (W) x 7.1" (T) x 6.1" (D)
 - 2-year warranty
 - 120V power supply

Item	Description	Price
SE03	Evacuator with filter, remote and adaptor, 120V	\$1,219.92
SEF06	Smoke Shark III Filter 2-Pack (non-sterile), 6 hr smart filter	\$117.63
SE95	SharkSkin smoke evacuation pencil adapter with 3/8" x 10' Tubing (20/bx, packaged sterile individually)	\$308.97



Bovie® Bantam PRO

Item: A952
Price: \$2,705.36
Max Wattage: 50
Bipolar: Yes
Cut: Yes
Voltage: 110-240
Electrode 3/32": Yes
 7.4" H x 8.98" W x 4.13" D
Weight (lbs): 5
Warranty (year): 4

Includes Accessories:
 Handpiece (autoclavable), Sheaths (*Non-Sterile x2, Sterile x2*), Grounding Cord & Pad, Wall Mount Kit, Electrodes - Sterile Disposable (*Blade x3, Needle x1, Ball x1, Blunt x2, Sharp x2, Derm-Elite™ x2*), Electrodes - Non Sterile Disposable (*Blunt x5, Sharp x5*), Power Cord & Manual.

Bovie® Derm 102

Item: A102
Price: \$744.38
Max Wattage: 10
Bipolar: Yes
Voltage: 110-240
Electrode 3/32": Yes
 7.2" H x 8.2" W x 2.7" D
Weight (lbs): 3
Warranty (year): 4



Includes Accessories:
 Handpiece (autoclavable), Sheaths (*Non-Sterile x2, Sterile x2*), Wall Mount Kit, Electrodes - Sterile Disposable (*Blunt x2, Sharp x2, Derm-Elite™ x2*), Electrodes - Non Sterile Disposable (*Blunt x5, Sharp x5, Derm-Elite™ x2*), Power Cord & Manual.

✓ Conmed Hyfrecator

Item: 7-900
Price: \$1,420.00
Max Wattage: 35
Bipolar: Yes
Voltage: 110
Electrode 3/32": Yes
 8.75" W x 7.5" H x 4" D
Weight (lbs): 6
Warranty (year): 1



Includes Accessories:
 Handpiece (autoclavable), Wall Mount Kit, Electrodes - Non Sterile Disposable (*Blunt x5, Sharp x5*), Power Cord & Manual.

✓ Bovie® Derm 942

Item: A942
Price: \$1,271.80
Max Wattage: 40
Bipolar: Yes
Voltage: 110-240
Electrode 3/32": Yes
 7.4" H x 8.98" W x 4.13" D
Weight (lbs): 5
Warranty (year): 4



Includes Accessories:
 Handpiece (autoclavable), Sheaths (*Non-Sterile x2, Sterile x2*), Wall Mount Kit, Electrodes - Sterile Disposable (*Blunt x2, Sharp x2, Derm-Elite™ x2*), Electrodes - Non Sterile Disposable (*Blunt x5, Sharp x5, Derm-Elite™ x2*), Power Cord & Manual.

For Bovie & Hyfrecator® Accessories go to delasco.com

Electrosurgery – Forceps




L to R: A827BP, A826, A824 & A825

Bovie® Bipolar Forceps

All forceps require a cord and foot switch ordered separately.

Item	Description	Price
A824	4½" Curved, Microfine, 0.5 mm	\$246.33
A825	4¾" Straight, Medium, 2 mm	\$246.33
A826	4½" Straight, Microfine, 0.5 mm	\$246.33
A827BP	Cord for A900/A940/A950	\$54.46
A827V	Cord for A102, A942, A952, 2250, 1250S & 3250	\$54.46



Delasco Bipolar Forceps

Top quality bipolar forceps and cord at a phenomenal price. Ideal for coagulation procedures in which stray electrosurgical currents are undesirable.

- Fits all standard electrosurgical units
- Forceps and cords are autoclavable


All forceps require a cord and foot switch ordered separately.

Item	Description	Price
26401	3½" Straight, McPherson, 0.5 mm	\$125.48
26402	4¾" Straight, Adson, 1.0 mm	\$125.48
26403	5½" Straight, Semkin, 0.5 mm	\$125.48
26404	6½" Straight, Cushing, 1.0 mm	\$125.48
26405	7" Straight, Cushing, 1.0 mm	\$125.48
26406	7" Straight, Cushing, 2.0 mm	\$125.48
26407	7" Straight, Gerald, 1.0 mm	\$125.48
26408	7¼" Bayonet, Gerald, 1.0 mm	\$125.48
26400	Standard Cord for Bovie, Conmed, and Most Other Electrosurgery Units	\$79.55

Fine Epilation Disposable Needles

- 18/8 stainless steel
- ⅜" Shaft


Item	Description	Price
DBEN	Sterile 0.004", 30/box	\$67.08



Bovie® Disposable Electrodes

- 3/32" Shaft

Item	Description	Price
ES02	Sterile Needle, 25/box	\$62.17
ES01	Sterile Blade, 50/box	\$78.94
ES22	Sterile 1/4" Loop, 5/box	\$71.20
ES24	Sterile 3/8" Loop, 5/box	\$71.20
ESSP	Electrode Scratch Pads, 40/box	\$67.28



Bovie® Reusable Electrodes

- 3/32" Shaft
- Autoclavable

Item	Description	Price
A830	Angled, Sharp, General	\$29.12
A831	Short, Angled, Ball	\$38.82
A834	Angled, Fine, Needle	\$42.06
ES22R	1/4" Loop, Tungsten	\$27.17



ConMed® Reusable Electrodes

- 3/32" Shaft
- Autoclavable

Item	Description	Price
711B	Angled Short Desiccation	\$40.92
714-S	Insulated Epilation Needle, 6/box	\$204.88

ConMed® Electrolase® Disposable Electrodes

Item	Description	Price
7-100-12	Non-sterile Sharp, 100/box	\$92.20
7-101-12	Non-sterile Blunt, 100/box	\$92.20
7-100-8	Sterile Sharp, 50/box	\$82.84
7-101-8	Sterile Blunt, 50/box	\$82.84



Delasco ElectroTips® Disposable Electrodes

Delasco ElectroTips are bent at the optimum operating angle and the hexagonal collar helps prevent twisting.

- 3/32" Shaft
- Color-coded insulation for easy identification

Item	Description	Price
✓ ET-101	Non-Sterile Sharp, 100/box	\$71.77
ET-102	Non-Sterile Blunt, 100/box	\$71.77
ET-103	Sterile Sharp, 50/box	\$85.28
ET-104	Sterile Blunt, 50/box	\$85.28



Delasco Electrolite® Disposable Electrodes

Physicians know that machined, spherical-radius tip electrodes are the preferred electrosurgical tool for precision coagulation with minimal tissue damage. The patented Electrolite disposable performs far above disposable flat-blade electrodes by giving you machined-quality technology in a disposable electrode.

- 3/32" Shaft
- Hexagonal collar prevents twisting of the electrode
- 3D shape typically found only in reusables
- Less unwanted tissue cutting
- Less concentrated debris buildup

Delasco Electrolite Guarantee: Try a box of Delasco Electrolite electrodes. If you aren't completely satisfied, call Delasco for full credit.

Item	Description	Price
✓ 26200	Non-Sterile, 100/box	\$81.06
26200S	Sterile, 50/box	\$73.17



Segal® Reusable Electrodes

- Made in the USA
- 3/32" Shaft
- Autoclavable

Item	Description	Price
✓ SE/L	Original	\$11.21
SE/E	Epilation Needle, 6/box (similar to the ConMed® 714)	\$104.19

Electrosurgery – Mobile Stand, Sheaths & Cautery




Delasco Universal Mobile Stand

Mounting board measures 6" x 7½". Pole is 37½" high. Base is 18½" wide, 4" high, and includes 5 swivel casters, 2 of which are locking.

- Delasco 734, 1134 and Electricator
- Bovie® units, all models
- Hyfrecators, all models
- Fits TCU wall mount bracket

Item	Description	Price
26800	Mobile Stand	\$253.08
26800/TC-PK	Tray & Clamp for Stand (Tray measures 13½" x 9¾")	\$166.94



Delasco ElectroSurgical Sheaths

A non-sterile sheath designed to fit over most electroSurgical handpieces. The use of handpiece sheaths significantly increases the life of the reusable handpiece and is a cost effective alternative to continuous sterilization.

We are so confident in this product, try them, if the quality and price of these don't exceed your expectations, send them back within 30 days for a full refund.

Not recommended for Hyfrecator 7-905 or 7-906 handpiece.

Item	Description	Price
✓ 26100	Non-Sterile 8", Latex-free, 100/box	\$23.53

Bovie® Handpiece Sheaths

(Compatible with Bovie and Hyfrecator handpieces.)

Item	Description	Price
A910ST	Sterile 38", Latex-free, 25/box	\$73.72

ConMed Protective Sheaths

(Compatible with Bovie and Hyfrecator handpieces.)

Item	Description	Price
7-796-18	Non-Sterile 8", Latex-free, 100/box	\$25.90
7-796-19	Sterile 36", Latex-free, 25/box	\$103.41

Disposable Cautery

Includes handle, tip, cover and non-replaceable batteries.

- Tips are not removable
- 2200° F
- Sterile

Item	Description	Price
✓ AA01	Fine Tip	\$12.52
AA03	Loop Tip	\$12.52

Change-a-Tip Cautery

Each tip includes one sterile handpiece sheath. Handle includes 2 replaceable AA batteries and 1 H101 tip.

- Sterile tips are interchangeable
- Tips have 4-year shelf life
- 2200° F

Item	Description	Price
HIT1	Power Handle w/H101 tip	\$56.25
H101	Fine Tip, 10/box	\$58.52
H103	Loop Tip, 10/box	\$58.52
H105	2" Flexitip (bendable tip), 10/box	\$66.65
ACS	Sterile Sheath 7½" L x 1¼" W, 10/box	\$30.01



MVE® Dewars

The MVE Lab Series of cryogenic dewars are high-efficiency, super insulated dewars that are the most convenient, economical way to store and dispense liquid nitrogen.

**Static Evaporation Rate is nominal. Actual rate may be affected by: nature of contents, atmospheric conditions, container history, or manufacturing tolerances.*

- Easy-to-maintain lid design
- High strength neck tube reduces liquid nitrogen loss
- Advanced Chemical Vacuum Retention System
- Advanced insulation system provides maximum thermal performance
- Lightweight with superior strength
- 5-year vacuum warranty

Item	MVE4	MVE5	MVE10	<input checked="" type="checkbox"/> MVE20	MVE30	MVE50
Price	\$988.98	\$893.58	\$942.34	\$1,222.18	\$1,449.02	\$2,131.38
Size (Liters)	4	5	10	20	32	50
SER* Liters/Day	0.19	0.15	0.18	0.18	0.22	0.49
Outside Diameter	7.3"	8.8"	10.3"	14.5"	17"	17"
Overall Height	16.8"	18.2"	21.5"	24.5"	24"	30.5"
Throat Size	1.4"	2.2"	2.2"	2.2"	2.5"	2.5"
Empty Weight (lbs)	6	8	12	19	25	31
Weight Full (lbs)	13	17	31	55	82	120



MVE® Manual Withdrawal Devices

- Dispenses 1-1½ cups of liquid nitrogen before losing pressure
- Allow for 10 minutes to rebuild pressure between uses
- May be left in place between uses



Brymill Cryo Tweezers

Designed for the localized treatment of skin tags.

- Extra bulk at the ends allows for longer freeze duration
- Narrow neck prevents cold from transferring down main body part of tweezer
- Made of stainless steel, brass & Teflon

Includes a Delrin Cup Holder and 10 Foam Cups.

Item	Description	Price
WD/MVE10	MVE10	\$971.61
<input checked="" type="checkbox"/> WD/MVE20	MVE20 (after 2007)	\$970.95
WD/MVE30	MVE30	\$970.95
WD/MVE50	MVE50	\$970.95

Item	Description	Price
BCT	Cryo Tweezer	\$293.33

Cryosurgery – Units & Accessories



CRY-AC® Straight Spray Needle

Item	Size	Price
CG/16/1	16G x 1"	\$36.00
CG/18/1	18G x 1"	\$36.00
<input checked="" type="checkbox"/> CG/20/1	20G x 1"	\$36.00



CRY-AC® Specialty Tips

Item	Size	Price
CG/AS	Acne Aperture	\$45.00
CG/T/E	Bullet Shaped E Tip w/Back Vent	\$31.32

CRY-AC® Replacement Tips

Item	Size	Price
CG/T/A	0.04"/1 mm	\$31.32
CG/T/B	0.031"/0.8 mm	\$31.32
CG/T/C	0.022"/0.55 mm	\$31.32
CG/T/D	0.016"/0.4 mm	\$31.32
CG/TS	Set of 4 (A, B, C & D)	\$112.28




CryoPro®

Durable precision - available in two sizes

350 ml

CryoPro® MINI
12-14 hrs. static holding time

500 ml

CryoPro® MAXI
20-24 hrs. static holding time

Buy CryoPro® for liquid nitrogen freezing treatment.

Cryosurgery is a minimally invasive treatment option that has proven both safe and effective within the field of dermatology for the treatment of various skin lesions such as actinic keratosis, viral warts, and basal cell carcinomas.

As cool as it gets



CryoPro® Liquid Nitrogen Apparatus

Offers full compatibility with other LN units which allows spray tips and contact probes from other systems to be used with both units.

Includes unit, 5 spray tips: A (1 mm), B (0.75 mm), C (0.55 mm), D (0.45 mm), 1 bent extension spray tip, 1 cleaning adapter and probe stand.

- Easy rotation for left or right hand users
- 3 year Warranty
- CE-marked



Brymill® CRY-AC®

Hand-held devices, these units set the industry gold standard in safety, versatility and control.

Includes unit, 4 spray tips: A (1 mm), B (0.8 mm), C (0.55 mm), D (0.4 mm) tips and one each 20 gauge x 1" straight and 20 gauge x 3" bent extension sprays.

- Easy rotation for left or right hand users
- 3 year Warranty
- CE-marked

Item	Size	Price
✓ CRYO/MAXI	16 oz, Static holding time 20-24 hours	\$831.76
✓ CRYO/MINI	12 oz, Static holding time 12-14 hours	\$831.76

Item	Size	Price
✓ CRYAC	16 oz, Static holding time 20-24 hours	\$810.54
✓ CRYAC3	10 oz, Static holding time 10-12 hours	\$810.54
CRYBABY	5 oz, Static holding time 3-4 hours	\$828.11



CryoOmega®

Clinically effective direct application of liquid nitrous oxide (-89°C/-127°F) in a convenient and economical disposable device. Cleared for 23 common lesions including verrucae, keratoses, and sebaceous hyperplasia. Most lesions require a single treatment. The internal cartridge provides 240 seconds (4 minutes) of total treatment time per unit.

- Convenient, hassle-free activation
- Can be used for multiple procedures

Item	Size	Price
3400	Single	\$445.05



Premier® Nitrospray Plus

A lightweight, self-pressurized liquid nitrogen device with a closed system that is trigger activated.

Includes unit, 5 spray tips (16, 17, 18, 19 and 20 gauge), a tip protector and base holder.

- Easy rotation for left or right hand users
- Luer Lock
- 3-year warranty
- Static holding time 8-12 hours

Item	Size	Price
6060	16oz	\$723.11
6065	10oz	\$723.11



Brymill
CRYOGENIC SYSTEMS
Established 1966

Freeze with ease!



Cryosurgery made easy.



Delasco
Preferred
Program

Ready to save? Join now.
delasco.com/preferred



Dynarex Sterile Towel Drapes

- 3-ply (Tissue/Poly/Tissue)
- 18" x 26", White

Item	Description	Price
DYN4409	Fenestrated, 50/box	\$19.00
DYN4410	Non-Fenestrated, 50/box	\$19.00



TIDI Table Goods

- Table Paper available in Choice or Everyday
- 99.5% protection against bacteria

Pillowcases

Item	Description	Price
919355	21" x 30", Fabricel®, 100/case	\$84.60

Drape Sheets

Item	Description	Price
918302	40" x 48", 100/case	\$38.30
918303	40" x 60", 100/case	\$48.90



Table Paper

Item	Description	Price
913182	"Choice" Smooth, 18" x 225', 12/case	\$70.45
913212	"Choice" Smooth, 21" x 225', 12/case	\$83.35
916183	"Choice" Crepe, 18" x 125', 12/case	\$57.12
913184	"Everyday" Smooth, 18" x 225', 12/case	\$65.11
913214	"Everyday" Smooth, 21" x 225', 12/case	\$78.06
916186	"Everyday" Crepe, 18" x 125', 12/case	\$44.93
916216	"Everyday" Crepe, 21" x 125', 12/case	\$52.59



TIDI Poly Tissue Capes and Gowns

- Tear-resistant, leak proof

Item	Description	Price
910415	30" x 21" 3ply, cape, white, 100/case	\$51.72
910420	30" x 42" 2ply, gown, white, 50/case	\$51.95
910520	30" x 42" 2ply, gown, blue, 50/case	\$51.70
910320	30" x 42" 3ply, white, 50/case	\$43.78
918507	30" x 42" 3ply, blue, 50/case	\$45.88

Professional Towel

- 13" x 18"
- White
- Poly/Tissue



Item	Description	Price
ALP1001A	3ply, 500/case	\$31.89
ALP1051	2ply, 500/case	\$36.24

Headbands

- Non-sterile



Item	Description	Price
GRA51823	2.5" W, White, Velcro, 48/pack	\$24.91
GRA82025	2.5" W, Black, Velcro, 48/pack	\$20.17



FlexFit® Exam Cover-Up

- Open sides, elastic band
- One size fits most

Item	Description	Price
GRA37263	25/case	\$67.11



NEW

Mayo Stand on Wheels

Includes Stainless Steel Tray 12-1/2" x 19".

- Plastic Mobile Base features locking casters
- Weight Capacity 11 lbs.
- Secure height adjustments 28" - 48"

Item	Description	Price
TEC4365	Stand and Stainless Steel Tray	\$191.74

- Cast aluminum 5 leg base
- Adjustable height from 29 1/2" to 49"
- Stainless steel 12 5/8" x 19 1/8" stationary tray
- Stainless steel 13 1/2" x 9 3/4" bottom adjustable tray

Item	Description	Price
<input checked="" type="checkbox"/> MS/WS	Stand with 1 Tray	\$249.08
MS/WD	Stand with 1 Tray & 1 Adjustable Tray	\$362.20
43431	MS/WD Adjustable Replacement Tray	\$168.44
43433	Replacement Clamp	\$117.27



NEW

Mayo Instrument Stand

- Stainless steel Heavy, weighted base.
- Telescoping rod adjusts tray to any height from 37 1/4" to 64"
- Includes Stainless Steel Tray 12 5/8" x 19 1/8"

Item	Description	Price
MS	Mayo Stand	\$444.52
MS/T	Stainless Steel Tray	\$90.58

Mayo Instrument Stand with U-Shaped Base

U-shaped base minimizes interference with other equipment

- Weight Capacity 11 lbs.
- Secure height adjustments 28" - 48"
- Includes Stainless Steel Tray 12 5/8" x 19 1/8"

Item	Description	Price
TEC4352	Stand and Stainless Steel Tray	\$178.62



Tutt-nauer T-EDGE11S T-Edge 11 Autoclave

The only gravity and airflow removal autoclave upgradable to pre/post vacuum. Designed to meet the most current sterilization standards ANSI/AAMI ST55 ensuring highest sterilization levels and highest quality construction to ensure longevity.

- More than 75% chamber capacity with 5 extra-long trays.
- Convenient and easy loading with door opening to 110°
- A variety of documentation and traceability options and barcodes.
- Displays programs, cycles times, water levels and more

Item	Description	Price
TEEDGE11	Tutt-nauer T-EDGE11S T-Edge 11 Autoclave	\$9,631.36

Tuttenauer T-Top 10 Autoclave

The T-Top10 autoclave is specifically developed for dental and medical clinics of small to medium practices, covering the fundamental needs for your general clinic sterilization with the aim of increasing the productivity of your practice. Designed to Meet Your Healthcare Clinics' Everyday Sterilization Needs.

T-Top offers universal capacity, handed to you, in an affordable and reliable solution you can count on for managing your daily practice sterilization process workload; an efficient autoclave which perfectly meets the requirements of all standard practices.

Efficiency with maximum reliability

- 21 Liter Large Chamber Capacity
- Class S Sterilizer (Closed-Door Drying)
- T-Connect: Advanced Remote Diagnostics
- Short Cycles

Item	Description	Price
TTOP10	21 Liter Chamber	\$6,364.00



Tuttnauer Valueklave 1730 Steam Sterilizer

Compact size makes it an ideal primary autoclave for small offices or secondary unit for larger practices.

- 2-gallon capacity
- Double lock to prevent door opening during cycle
- Dual thermostat to protect against overheating
- Automatic shut-off
- Unwrapped cycle time: cold, 12 minutes; hot, 9 minutes
- International certifications

Item	Description	Price
V1730	Tuttnauer Valueklave 1730 Steam Sterilizer	\$3,488.10



Tuttnauer's Clean & Simple Ultrasonic Cleaner with Basket

Provides 60 kHz of quiet cleaning power which is less aggressive to instruments while still removing all surface contamination.

- Stainless steel tank with cover
- 1-Gallon Tank Dimensions 9 3/8" W x 4" H x 5 3/8" D
- 1-Gallon External Dimensions 10 1/4" W x 8 1/4" H x 6 1/2" D
- 3-Gallon Tank Dimensions 15 1/2" W x 6.0" H x 9.0" D
- 3-Gallon External Dimensions 16 1/2" W x 12 1/8" H x 10.0" D
- Includes stainless steel self draining basket

The 1-gallon model holds a standard selection of instruments or 2 half sized cassettes. The 3-gallon model holds 3 full size cassettes.

Item	Description	Price
CSU1BK	1-Gallon, 110 Volt	\$864.80
CSU1HBK	1-Gallon, Heated, 110 Volt	\$978.33
CSU3BK	3-Gallon, 110 Volt	\$2,020.38




Tuttnauer Cleaners

- Convenient, pre-measured packets
- Cleans chamber, tubing and reservoir
- Uses less storage space
- Replaces messy liquid cleaners

Item	Description	Price
CB0010	10 packets/box	\$70.13
CS0144	144 tablets/box	\$126.21



Wound Care – Bandages, Skin Closure & Irrigation



Leukoplast Fabric Bandages

- Fabric
- Sterile
- Long-lasting adhesive
- Latex-free

Item	Description	Price
7645611	4-Wing 3" x 3", 50/box	\$19.55
7645616	Spots, 7/8", 100/Box	\$8.05
7645610	Patches 2" x 3", 50/box	\$11.50
7645612	Patches, 1 1/2" x 2", 100/Box	\$13.23
7645606	Strips 1" x 3", 100/box	\$9.20
7645617	Strips 3/4" x 3", 100/box	\$8.05



3M™ Steri-Strip™

- Porous, non-woven elastic backing
- Pressure-sensitive, hypoallergenic adhesive skin closures
- White strips include polyester filaments for added strength
- Blend tone is non-reinforced allowing natural skin color to show
- Sterile

Each box includes 50 envelopes of 3, 5 or 6 strips.

Item	Description	Price
✓ R1540	White, 1/8" x 3", 5 strips/envelope	\$83.11
R1541	White, 1/4" x 3", 3 strips/envelope	\$83.11
R1542	White, 1/4" x 1 1/2", 6 strips/envelope	\$83.11
1551	Blend Tone, 1/4" x 3", 3 strips/envelope	\$85.31
B1557	Blend Tone, 1/2" x 4", 6 strips/envelope	\$142.16



American White Cross Fabric Bandages

- Highly absorbant non-stick pad
- Long-lasting adhesive
- Sterile

Item	Description	Price
SVT3D	Strips 3/4" x 3", 100/box	\$4.48
SV13D	Strips 1" x 3", 100/box	\$4.48
OVKXD	Spots 7/8" round, 100/box	\$4.48
SVWXD	Spots 1 1/4" oval, 100/box	\$4.48
PVJ2D	Patches 1 1/2" x 2", 100/box	\$6.73
PV23R	Patches 2" x 3", 50/box	\$5.60



Dukal Wound Closing Strips

- Non-woven reinforced fibers
- Translucent and porous
- Sterile and latex-free

Each box includes 50 envelopes of 4, 5, 6, or 10 strips.

Item	Description	Price
✓ DUK5150	White, 1/8" x 3", 5 strips/envelope	\$44.82
DUK5156	White, 1/4" x 4", 10 strips/envelope	\$63.86
DUK5157	White, 1/2" x 4", 6 strips/envelope	\$63.86
DUK5158	White, 1" x 4", 4 strips/envelope	\$63.86

SAVE THE DATE



**AMERICAN
DERMOSCOPY
MEETING**

JUNE 26-28, 2025
BONUS WORKSHOP JUNE 25TH

MONTEREY MARRIOTT
MONTEREY, CA

NEAR PINNACLES NATIONAL PARK

REGISTER AT AMERICANDERMOSCOPY.COM

Mastisol® Adhesive

- Liquid adhesive
- Clear, non-irritating and non-water soluble
- Secures the most difficult dressings over a prolonged period of time
- Latex-free

Item	Description	Price
<input checked="" type="checkbox"/> MAST/1	2/3 ml Applicator	\$4.59
MAST/1/2	½ oz	\$22.95
MAST/2	2 oz Dispenser Cap	\$57.38 ◆



Wound Care – Sponges, Wraps & Dressings



Dukal Woven Gauze Sponges

Ideal for general cleaning, dressings, prepping and debriding wounds.

- Made from 100% premium bleached woven cotton
- Highly absorbent

Item	Description	Price
DUK6412	4" x 4", 12 ply, 2/pack, 25 pk/bx, Sterile	\$4.56
DUK416-10	4" x 4", 16 ply, 10/tray, Sterile	\$1.77
DUK8505	2" x 2", 8 ply, 200/sleeve, Non-sterile	\$1.12
DUK8506	2" x 2", 12 ply, 200/sleeve, Non-sterile	\$1.72
DUK8508	3" x 3", 12ply, 200/sleeve, Non-sterile	\$3.73
DUK8509	4" x 4", 8 ply, 200/sleeve, Non-sterile	\$4.87
DUK8510	4" x 4", 12 ply, 200/sleeve, Non-sterile	\$6.83
DUK4162	4" x 4", 16 ply, 200/sleeve, Non-sterile	\$11.89



Ethicon Surgiform®

Used dry or saturated with sterile sodium chloride solution, is indicated for surgical procedures (except ophthalmic) for hemostasis, when control of capillary, venous and arteriolar bleeding by pressure, ligation and other conventional procedures is ineffective or impractical.

- Sterile
- Water-insoluble, malleable, porcine gelatin absorbable sponge
- Hemostasis agent

Item	Description	Price
ETH1972	12 sq cm x 7 mm, 12/box	\$205.33



Dukal Non-Woven Sponges

Highly absorbent Rayon/Poly blend, Virtually non-linting multi-purpose sponges. Ideal for general wound care and applying ointments.

- Absorbent, non-linting and soft multi-purpose sponge
- Rayon/poly blend

Item	Description	Price
DUK2102	2" x 2", 4 ply, 200/sleeve, Non-sterile	\$1.16
DUK2103	3" x 3", 4 ply, 200/sleeve, Non-sterile	\$2.45
DUK2104	4" x 4", 4 ply, 200/sleeve, Non-sterile	\$4.32



Hypafix® Dressing Retention Sheet

- Non-woven, polyester fabric
- Layer of non-irritating adhesive
- Secure dressings, tubes, catheters and instruments
- Stretchable
- Hypoallergenic
- Water-resistant

Item	Description	Price
25211	2" x 2 yards	\$8.41
25212	2" x 11 yards	\$16.80
25213	4" x 11 yards	\$25.21
25214	6" x 11 yards	\$38.28



Dukal Cohesive Elastic Bandage

- Non-sterile
- Self-adherent elastic bandage
- Contains Latex

Item	Description	Price
DUK8025T	2" x 5 yards, 36 rolls/box	\$65.29
DUK8035T	3" x 5 yards, 24 rolls/box	\$52.50
DUK8045T	4" x 5 yards, 18 rolls/box	\$50.22

Burton™ UV LED Magnifier Lamp

Used for inspection for detecting defects. A lightweight alternative to other Woods lamps. Can be used with cord or go cordless and rely on battery-saving sleep mode.

- 1.85" Height x 3.85" width viewing area
- 3,288 Lux at 6" (White LED)
- 95 CRI, 3800 Kelvin (White LED)
- 36,000+ hours (White LED)
- 30,000+ hours (UV LED 365nm)
- Easy to clean and disinfect
- 4 hours of continuous use

Item	Description	Price
31604	100-240V rechargeable	\$626.24



Epic LED Magnifier

LED magnifier with light quality that creates a large viewing area with an arm and neck system that is rugged and flexible.

- 1.88x magnification (3.5 diopter)
- Large optical viewing area
- Superb light quality; 90 CRI, 4000 Kelvin
- 10-100% dimming capability
- Auto shut-off
- Low maintenance LED array
- Medical standards approved
- 5-year limited warranty

Item	Description	Price
EMLFL	Floor Model	\$997.00
EMLW	Wall Model	\$760.00



Burton Nova Exam LED Light

- The LED source provides a virtually heat-free light so the patient remains comfortable.
- Dual handles are easily accessible to position the light head.
- The patented K-arm moves smoothly with no drifting.
- Maximum arm reach of 45" (114.30 cm)

Item	Description	Price
NXFL	Floor Stand	\$816.00



Burton Medical UV501 Woods Lamp (no magnifier)

- Special fluorescent bulbs produce UV-A at approximately 360nm wavelength, not harmful to eyes or skin.
- Includes: 2 UV bulbs
- 6,000 average hour bulb life
- Assembled in the USA
- 5-year manufacturer warranty

Item	Description	Price
31501	handheld, 2 UV bulbs	\$438.99



Burton Nova Exam LED Light

- The LED source provides a virtually heat-free light so the patient remains comfortable.
- Dual handles are easily accessible to position the light head.
- The patented K-arm moves smoothly with no drifting.
- Maximum arm reach of 45" (114.30 cm)

Item	Description	Price
NXC	Celing Mount	\$816.00



Burton Outpatient (Wall Mount)

- 87, lux (8 footcandles) at 24" (61cm)
- 33 K color temperature
- CRI (Color Rendering Index) of 94
- Focuses with central removable, autoclavable SteriHandle™
- Friction knob allows user to adjust friction on the self-balancing arm positioning to desired tension without a tool
- Cool operation assured with heat filters and dichroic coatings

Also available in wall mount, double ceiling and 220 volt by special order.

Item	Description	Price
224142	Wall Mount	\$3,245.00



Burton Medical Super Exam Light, LED (Floor Stand)

- 24,000 Lux at 39.37" (1 m)
- CRI of 95
- Adjustable color temperature of 3500K/4000K/4500K
- Lower color temperature emphasizes the reds, while the higher temperature provides cooler white for skin tones
- User-friendly control panel for color temperature
- Dimming feature 100% - 0%
- Intuitive and ergonomic for one-handed operation

Item	Description	Price
SELEDACFL	Floor Stand	\$1,319.00



Burton Outpatient LED Exam Light 110 Volt (Floor Stand)

- 90,000 lux (8361 fc) at 24" (61 cm)
- Color Temperature: 3700 K
- CRI (Color Rendering Index): 95
- Light Emitting Diode (LED)
- 20,000 hours average LED life
- Power: 42 Watts
- Power Supply: Integral electronic driver
- Mounting Options: Floorstand, Wall, Ceiling
- Maximum Arm Reach: 37.5" (95.25 cm)

Item	Description	Price
OP316FL	Floor Stand	\$3,436.00



Burton AIM LED Exam Light Universal Input (Floor Stand)

- 45,000 lux (4181 fc) at 1 meter
- 4300 K color temperature
- CRI (Color Rendering Index) of 92
- 3 LED modules with a 20,000-hour life
- Adjustable focus using center handle
- UL60601-1 / IEC 60601-1 / IEC 60601-1-2 / IEC 60601-2-41; CAN/CSA C22.2 601-1 M90 certified
- Assembled in the USA
- Universal Input 100-240V
- Standard US Plug, Nema 5-15P

Item	Description	Price
ALED/FL	Floor Stand	\$4,158.00
ALED/SC	Ceiling Mount	\$4,222.66

L'Applique™ Lotion Applicator

Allows patients to apply medications, creams, lotions, washes or moisturizers without stretching or twisting.

- Applicator and pads are latex-free
- Includes one attached pad

Item	Description	Price
✓ LA	Ivory	\$7.96
LA/RD	Red	\$8.28
LA/PR	Purple	\$8.28
LA/P	Replacement Pads, Pair	\$2.56



Gel Packs

Gel masks relieve post-surgery pain and swelling. Commonly used for dermal fillers, facial injections, chemical peels, blepharoplasty, oral surgery and other procedures.

- Remains flexible while frozen
- Freezes quickly for reuse
- Cloth surface provides a comfortable protective barrier
- Non-toxic
- Cannot be warmed or heated

Item	Description	Price
GEL/RECT	Rectangle, 4" x 10"	\$5.12
✓ GEL/ROUND	Round, 4"	\$4.56



Dynarex Instant Cold Pack

Single squeeze allows for easy activation of packs, offering immediate cooling.

- No refrigeration needed
- Flexible material conforms to body contours
- Non-Toxic

Item	Size	Price
DYN4517	4" x 5", 24/case	\$16.62



Aquaphilic Ointment

Hydrated hydrophilic ointment. A superior oil in water lubricating cream base, free from fragrance, including masking. Originally formulated at the University of Iowa for the asteatosis associated with dry mid-western winters.

Item	Size	Price
✓ AQUA	1lb	\$16.53
AQUA/U	10% Urea, 1lb	\$16.71
AQUA/U/20	20% Urea, 1lb	\$16.71



Ointment Jars


Pack of 10 white jars in various sizes.
Made of polypropylene plastic.

Item	Size	Price
OJ/0-10	½oz	\$12.77
OJ/1-10	1oz	\$8.40
OJ/2-10	2oz	\$20.17
OJ/4-10	4oz	\$11.21
OJ/16-10	16oz	\$20.17



Cotton Applicators

Item	Size	Price
APP/C/6	6" wood 1000/pack	\$11.44
AMS9841	Puritan, 6" wood Sterile, 50/pack	\$7.53
806WCL	6" wood 500/pack	\$16.38



Disposable Fan Brushes

Item	Size	Price
511013PK	7¾"L, 1½" W Taklon, 4/pack	\$10.08



Hand-Held Mirror

Item	Size	Price
651BK	9" x 8"	\$6.62



Dermal Gloves

- Lightweight 100% cotton
- Provides stain protection from hand creams, ointments, and medication

Item	Size	Price
9665	Regular (Small/Medium), 12 pair/pack	\$25.04
9666	Large (Medium/Large), 12 pair/pack	\$26.92



Clear Jar

5 gram jar. Clear jar with threaded cap.
Made of polystyrene.

Item	Size	Price
CJ5	5gm/100 pack	\$27.28

CNH® Ear Protector

Useful to relieve pressure on the ear after ear surgery, or to relieve the discomfort associated with Chondrodermatitis Nodularis Helicis (CNH) also known as Painful Nodules of the Ear.

The CNH ear protector covers the ear and protects it from pressure while turning during sleep.

- Specially designed foam is flexible yet very dense
- May be washed in mild soap and air-dried
- Latex-free
- Velcro strap

Not intended to relieve sustained, direct pressure on the affected ear but meant for inadvertent pressure as when sleeping.



Item	Description	Price
✓ CNH-1	Ear Protector	\$22.73



Dr. Dan's CortiBalm™

- #1 dermatologist-recommended lip balm
- Effective on severely chapped lips
- Contains beeswax, petroleum jelly, mineral oil and 1% hydrocortisone
- Not for use on cold sores or fever blisters

Item	Size	Price
CBALM	.14oz Tube	\$7.23
JCBALM	.25oz Jar	\$7.23



ALO-AR® Essential Oil

- Oil used in treatment of alopecia areata
- Includes grapeseed, rosemary, lavender, thyme and cedarwood oils

Item	Size	Price
ALO	1oz	\$40.92

Cosmetics – Ear Piercing



Coren® PS Earpiercer

- Firm control
- Not autoclavable
- Sterile
- All earring types included with the Coren Ear Piercer do contain trace amounts of nickel

Item	Description	Price
COREN-SQ	Grey Flex Handle Only	\$5.63
COREN-SS	Piercer w/Stainless Steel Earrings, 6 pair	\$56.30
COREN-GOLD	Piercer w/24k Gold Plated Earrings, 6 pair	\$56.30
COREN-G/D	Piercer w/24k Gold Plated Simulated Diamond, 6 pair	\$58.54
COREN-PEDI	Piercer w/Childrens 24k Gold Plated Earrings, 6 pair	\$61.92
COREN-SS/D	Piercer w/Surgical Stainless Steel w/ Simulated Diamond, 6 pair	\$56.30

Blomdahl® Medical Ear Piercing Kit

Includes ear piercing gun, 10 medical plastic single earring cassettes, 10 antiseptic wipes, a surgical marking pen, hair clip, mirror, dummy ear lobes, promotional material, consent forms and registration forms.

- Uses sterile, disposable cassettes
- Instrument does not contact ear
- Earrings included are medical grade plastic, 4 mm diameter

Item	Description	Price
BLOM/SK	Starter Kit w/Accessories	\$321.08
BLOM/REG	Regular Ear Piercing Gun	\$107.80




Blomdahl® Mini Ear Piercing Gun

Includes the gun, one surgical marking pen, one medical plastic single earring cassette and one antiseptic wipe.

- Same size as regular ear piercing gun
- Utilizes smaller pin shaft for mini earring cartridges
- Titanium earrings are nickle-free, 3 mm diameter

Item	Description	Price
BLOM/MINI	Blomdahl Mini Kit w/Accessories	\$103.82
MINI/GTC	Golden Titanium, Crystal, 10/box	\$106.89
MINI/TBC	Silver Titanium, Crystal, 10/box	\$106.89






MediSystem Ear Piercing

- Fits regular and mini earrings, is easy to use and reliable
- Made of solid construction with hardened stainless steel working parts and injected molded handle
- Not autoclavable

Item	Description	Price
MS1000N	Ear Piercing Gun	\$76.61
MS2010	24K Gold Plated Ball, 12 pair	\$21.08
MS3010	Stainless Steel Ball, 12 pair	\$21.08
MS2090	24K Gold Plated Cubic Zirconia, Tiffany Mount, 12 pair	\$64.40
MS3090P	Stainless Steel Cubic Zirconia, Tiffany Mount, 12 pair	\$64.40
MS2010M	Mini 24K Gold Plated Ball, 12 pair	\$19.91



Pierced Ear Solution

- Anti-bacterial, germicidal topical cleanser for pierced ear area
- Promotes healing
- Contains 0.1% benzalkonium chloride

Item	Description	Price
MS1003	4oz Bottle, 24/pack	\$77.28

Blomdahl® Regular Earring Cassettes

- Titanium earrings are nickel-free and hypo-allergenic
- Silver titanium are rhodium plated
- Golden titanium are plated with 24 karat gold on surfaces not in contact with the skin
- Single earring per cassette

Item	Description	Price
BLOM/MP/41PG	Med. Plastic Crystal, Daisy, 10/box	\$106.89
BLOM/MP/42PG	Med. Plastic, Crystal Daisy Rose, 10/box	\$106.89
BLOM/MP/45PG	Med. Plastic Crystal, Daisy Rainbow, 10/box	\$106.89
BLOM/MP/46PG	Med. Plastic Crystal, Daisy Alex/Rose, 10/box	\$106.89
BLOM/GTC	Golden Titanium Crystal, 10/box	\$106.89
BLOM/GTP	Golden Titanium Plain, 10/box	\$106.89
BLOM/TBC	Silver Titanium Crystal, 10/box	\$106.89
BLOM/GTC/T	Golden Titanium Crystal Tiffany, 10/box	\$134.69
BLOM/TBC/T	Silver Titanium Crystal Tiffany, 10/box	\$134.68
BLOM/PRL	Golden Titanium Pearl, 10/box	\$134.68
BLOM/STP	Silver Titanium Pearl, 10/box	\$134.68
BLOM/MP	Medical Plastic, 10/box	\$91.17

Protection – Gloves & Disinfectants

Nitrile Non-Sterile Exam Gloves



Ventyv® Elephant

Textured, 3.5 mil, Violet Blue Color.

- Chemo approved
- 200/box

Price \$9.24

10333101	XS	10335101	M
10334101	S	10336101	L



Ventyv® Walrus

Finger Texture, 3 mil, Cobalt Blue

- Chemo approved
- 300/box (XL - 200/box)

Price \$13.24

10346102	S	10348102	L
10347102	M	10349102	XL



Ventyv® Stallion

Finger Texture, 3.5 mil, Black.

- Chemo approved
- Fentanyl Tested
- 100/box

Price \$8.90

10319107	S	10321107	L
10320107	M		



Ventyv® Flamingo

Finger Texture, 3 mil, Coral.

- Formulated with Aloe & Vitamin E
- 200/box (XL - 180/box)

Price \$10.00

10333106	XS	10336106	L
10334106	S	10337106	XL
10335106	M		



Ventyv® Kangaroo

Textured, 4 mil, Blue.

- Chemo approved
- 200/box

Price \$10.86

10333103	XS	10335103	M
10334103	S	10336103	L



Ansell® Microtouch®

Finger Texture, 3.1 mil, Pink.

- Chemo approved
- 100/box

Price \$20.70

6034510	XS	6034512	M
6034511	S	6034513	L

Latex Non-Sterile Exam Gloves



Ventyv® Polymed

Finger Texture - 3.15 mil
 Natural White - 100/box
 • Polymer Coated, Powder-Free

Price	\$7.83			
PM101	XS	PM103	M	
PM102	S	PM104	L	

Latex Sterile Surgical Gloves



Ansell® Sensi-Touch

Smooth Texture - 6.7 mil White - 50 pair

Price	\$60.03			
7821PF	5.5	7824PF	7	
7822PF	6	7825PF	7.5	
7823PF	6.5	7826PF	8	



CaviWipes™ HP Disinfecting Wipes

The newest member of the trusted CaviWipes family, CaviWipes HP provide a fast 1-minute universal contact time against 67 pathogens including MDROs.

- Alcohol-free hydrogen peroxide formulation

Item	Description	Price
MET16-1100	160 Wipes, 6" x 6.75"	\$16.37
MET16-1150	65 Wipes, 9" x 12"	\$19.61

Neoprene Sterile Surgical Gloves



Ansell® Gammex

Smooth Texture - 7.5 mil Green - 50 pair

Price	\$196.65			
8512	6	8515	7.5	
8513	6.5	8516	8	
8514	7	8517	8.5	

CaviWipes™ Disinfecting Wipes & Surface Disinfectant

CaviWipes™ and CaviCide™ are intermediate level disinfectants recommended for use on hard, non-porous surfaces.

- All in 1 cleaner and disinfectant
- 3 minute kill time: (Both) TB, MRSA, (CaviCide™) VRE, S. Aureus, Pseudomonas aeruginosa, S. Enterica & T. Mentagrophytes
- Effective against (Both) HIV-1, HBV, and HCV, (CaviCide™) HSV-1, HSV-2, Influenza A2 Virus in 2 minutes

Item	Description	Price
MET13-1100	160 Wipes	\$16.50 G
MET13-1024	24 oz Spray	\$15.84 G
MET13-1000	Gallon	\$41.19 G





Durette II/III Patient Laser Shields

Includes 1 pair eyeshields, 2 elastic headbands and black case.

- 1 mm thick stainless steel cups
- Regular, 28 mm x 49 mm
- Autoclavable


Item	Description	Price
21-335	Durette II Patient Laser Shields	\$316.17
21-336	Durette II Replacement Headband	\$51.96
21-440	Durette III Nasal & Temporal Kit	\$111.00
21-445	Durette III Patient Laser Shields	\$309.01

Noir® Stainless Steel or Polymer IPL/LED/Laser Eye Shield

Includes case and lanyard.

- Stainless steel or fiber reinforced polymer
- For OD 7+ @ 190-10,600 nm
- Adjustable nosepiece
- CE approved

Item	Description	Price
I-SHIELD	Eye Shield	\$149.01
P-SHIELD	Eye Shield	\$69.68
I-SHIELD/R	Replacement Silicone Ring Pair	\$11.15
I-SHIELD/S	Replacement Strap	\$5.76



Ocular Protective Shields

- For conventional lid surgery
- Autoclavable with gas or steam
- Life of 50-100 autoclavings

Not for use with laser. Includes an autoclavable suction cup for easy removal. Cold sterilization is not recommended.

Item	Description	Price
PS/1P	X-Small, 23.5 x 21.5 mm, one shield	\$89.87
PS/2P	Small, 26 x 23.5 mm, one shield	\$89.87
PS/3P	Medium, 27.5 x 24.5 mm, one shield	\$89.87
PS/4P	Large, 28.5 x 25.5 mm, one shield	\$89.87
23-012LF	Suction Cup, 12/pack	\$48.88



Theia Eye Block Disposable Eye Shields

- Secure lock tab ensures best protection and prevents lifting
- Non-adhesive tab allows for easy application and removal
- Special formulated adhesive to secure in key areas and eliminate light leakage
- Hypoallergenic

Laser Eye Shield

- For use with Laser, IPL and LED

Comparable Honeywell Laser-Aid™ Disposable Eye Shield.

Item	Description	Price
<input checked="" type="checkbox"/> 31-8300	50/box (25 pair)	\$91.00

IPL Non-Laser Eye Shield

- For use with IPL and microdermabrasion
- IPL-Eye Patch Tested for IPL visible and IR range of 500-1,100nm

Comparable to Honeywell IPL-Aid™ & Derm-Aid™ Disposable Adhesive Eye Shield.

Item	Description	Price
<input checked="" type="checkbox"/> 31-9300	IPL Eye Shield, 100/box (50 pair)	\$95.50

Derm-Eye Shield

- for use with microdermabrasion and facials
- Not suitable for use with lasers, IPL or LED

Item	Description	Price
DA-100	Derm Eye Shield, 100/box (50 pair)	\$68.00

Delasco Disposable Eye Shields



- **FDA Cleared Plus Savings:** FDA Class 1 Medical Devices, offering effective protection at a lower price than other brands.
- **Comfort-Focused Design:** Made from soft, hypoallergenic materials, these shields ensure a snug fit without causing irritation.
- **Prevents Lifting:** The innovative lock tab secures the eye shield in place, preventing any lifting during use.
- **Easy Application and Removal:** Featuring a non-adhesive tab, these shields allow for quick and hassle-free application and removal.
- **Color Coded for Easy ID:** Each type of eye shield is marked with a unique color on the top and side of the box, minimizing the risk of selecting the wrong product.



Derm Eye Shield

- Protects eyes during standard dermatological procedures like microdermabrasion

Item	Description	Price
31-7300D	100/box (50 pair)	\$56.05



IPL Eye Shield

- Protects eyes during IPL and LED procedures
- Protection from IPL 420-1200 nm

Item	Description	Price
31-9300D	100/box (50 pair)	\$78.85



Laser Eye Shield

- Protects eyes during laser, IPL, and LED procedures
- Protection from medical lasers operating OD7 < 1064 nm and OD6 < 10,600 nm
- Protection from IPL 420-1200 nm

Item	Description	Price
31-8300D	50/box (25 pair)	\$75.05

UV Goggles

- Complete protection during UVA or UVB treatments
- Latex-free

Item	Description	Price
DPE-1	12/pack	\$90.75



Cox II Patient Laser Corneal Eye Shields

For use during laser, IPL and blepharoplasty. Includes a suction cup for easy removal.

- Shields and suction cups can be sterilized in solutions, gas or autoclaved
- Sold individually, not in pairs
- Stainless steel shields
- CE marked

Item	Description	Price
CS/5	XS, 22 x 19.5 mm, each	\$140.68
✓ CS/6	S, 24 x 21 mm, each	\$140.68
CS/8	M, 25.5 x 23 mm, each	\$140.68
CS/7	L, 27 x 24 mm, each	\$140.68





4001 E Plano Parkway
Suite 100
Plano, TX 75074-1823

25WSA

PRSR-T STD
US POSTAGE
PAID
DELASCO

All prices represent retail pricing and are subject to change. Items are subject to availability and may not be available at Delasco.com.

Save your practice
time & money. Join our
**Preferred
Program**

Use the QR code to find out more



Delasco Biopsy Punches

Explore surgery products starting
on page 41



**Delasco
Has
Everything
You Need**

See our full selection of
products at delasco.com

Delasco Chemicals
Find more chemicals
starting on page 05

